

DermaPlast® EFFECT Corn Plasters

General Product Description

DermaPlast EFFECT Corn Plasters presentations are composed of a layer of 300 µm / 0,3 mm hydrocolloid material acting as the wound pad, another larger layer of 300 µm / 0,3 mm hydrocolloid material acting as the dressing backing material, the silicone paper protection and the peel pack (individually sealed plaster). All individually sealed plasters (one plaster in each peel pack) are packed in a labeled plastic box, which includes a leaflet with instructions for use.

Intended Purpose

Hydrocolloid plasters are non-active, single-use, self-adhesive, non-sterile devices consisting of a self-adhesive hydrocolloid wound contact layer and an adhesive coated backing material (carrier material), protected with a silicone paper, for the protection and treatment of blisters, corns, minor injuries/cuts, cold sores, abrasions and absorption of exudates.

The products are intended for humans, without restrictions (in age group, weight range, health status or condition), for use by lay-users, in contact with intact skin (healthy skin) and breached skin, and short-term (60 minutes ≤ Duration ≤ 30 days, acc. to MDR) / prolonged use (24 hours < Duration ≤ 30 days, acc. to ISO 10993-1:2018) duration of use.

The products do not contain resorbable components, medicinal products / drug substances (derived from human blood or plasma or not), tissues or cells of human or animal origin, radioactive substances, nanomaterials or other hazardous components.

Application/Indication

- Wash hands and ensure the affected area is clean, dry and oil free before applying the plaster.
- Open the peel pack and extract the plaster.
- Warm the plaster with the hands before applying.
- Remove the first protective paper and apply the exposed part of the plaster to the affected skin area.
- Then gently pull the second protective paper and apply the rest of the plaster. Gently mould the plaster edges to ensure optimum adhesion.
- The application maximum of the hydrocolloid is up to 4 days per individual product. Should you wish to remove the plaster early, first wet it briefly with warm water.

Reference Numbers

Ref no.	Description	Length	Content of folding box
5365401	DermaPlast EFFECT Corn Plaster	40x17mm	P7

DermaPlast® EFFECT Corn Plasters

Residual Risks, Contra-Indications and Undesirable Side Effects, Warnings

The hydrocolloid plasters must not be used on non-treated infected corns, major injuries, highly bleeding or exudative wounds or for pressure ulcers.

All presentations:

- Only apply the plaster if the hydrocolloid pad completely covers the blister/minor injury/corn/cold sore/abrasion.
- In cases of diabetes mellitus always consult a doctor before using plasters.
- Do not apply the product for more than 4 days.
- Store plasters in a cool (not in the fridge), dry place, out of the reach of children.

All presentations except Cold sore presentations:

- Consult a doctor before applying the plaster on a potentially infected blister, minor injury, or if the blister or minor injury shows signs of infection while plasters are in place (e.g. redness, overheating, swelling or pain).

Specially for Blister presentations:

- Do not burst a sealed blister as this could risk an infection.

Single Use Device

Reusing a disposable and single-use medical device is dangerous and it could cause risk of infection. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Product Disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components of DermaPlast® EFFECT Corn Plasters should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Product Performance /Material Characteristics

Formulation Hydrocolloid 300 µm/0,3 mm			
Property	Test method	Unit	Specification
Peel adhesion on Steel Plates (90°)	TM00 0090-04	N/25mm	6,9 (+2,7 / -4,1)

DermaPlast® EFFECT Corn Plasters

Labelling

Lot-No. with 9-Digit Code

e.g.:



0XXX

Year and number
of production

XX

Week of
production

12X

For internal purposes
only

Manufacturing Date

e.g.:



2020
year

01
month

07
day

Use-by-Date

e.g.:



2023
year

-

01
month

Shelf Life: 3 years

Medical Device



Unique Device Identification (UDI)



Distributor



Distributor

Latest Date of Revision: 2020-03-12