RFslide

A RAPID SLIDE TEST FOR THE QUALITATIVE AND QUANTITATIVE DETERMINATION OF RHEUMATOID FAKTOR IN HUMAN-SERUM

Principle:

The RF reagent is a suspension of polystyrene latex particles sensitized with specially prepared human IgG in order to avoid inespecific agglutination . The RF latex reagent sensitivity has been adjusted to detect a minimum of 8IU/ml of rheumatoid factors according with the WHO international standard without previous sample dilution.

Summary:

Rheumatoid factors are a heterogeneous group of autoantibodies directed against the antigenic determinants on the Fc-region of IgG molecules. They are important in the diagnosis of rheumatoid arthritis, but can also be found in other inflammatory-rheumatic diseases and in various non-rheumatic diseases. They are also found in clinically healthy persons over 60 years of age. Despite these restrictions, the detection of rheumatoid factors is a diagnostic criterion of the American College of Rheumatology for classifying rheumatoid arthritis. These autoantibodies occur in all the immunoglobulin classes, although the usual analytical methods are limited to the detection of rheumatoid factors of the IgM type.

The classic procedure for the quantitation of rheumatoid factors is by agglutination with IgG-sensitized sheep erythrocytes or latex particles. Particular problems of these semi-quantitative methods are the poor between-laboratory precision and reproducibility, together with standardization difficulties. For these reasons, new assay methods such as nephelometry, turbidimetry, enzyme-immuno-assays and radioimmunoassays have been developed. The Biocon RF assay is based on the immunological agglutination principle with enhancement of the reaction by latex particles.

Reagent Concentration:

R1:

Latex particles coated with anti-RF preservative

R2:

Protein containing solution with rheumatoid factors > 8 IU/ml preservative

R3:

Protein containing solution with rheumatoid factors < 8 IU/ml preservative

Preparation and stability

The reagents and control sera are stable:

Up to the expiration date at +2°C to +8°C Do not freeze!

Specimen:

Collect serum using standard sampling tubes
Use fresh clear serum specimens. Plasma, lipemic serum or microbial contamination may cause erroneous results.
If the test cannot be performed immediately, store the specimen at +2 to +8°C for up to 72 h.
For longer storage, freeze the serum.
Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
Exercise the normal precautions required for handling all laboratory reagents.
Contaminated sera and a linger reaction time than 3 min. may cause false positive agglutination.

Precautions:

Reagents containing sodium acid may combine with copper and lead plumbing to form highly explosive metal azides. Dispose of reagent by flushing with large amounts of water to prevent azide build up. The positive and negative controls were prepared from human sera which have been tested using FDA licensed methods and found to be non reactive for HbsAg and HIV antibodies. However, no test method can offer complete assurance that infectious agents are absent. Therefore all human specimens should be considered potentially infectious.

Procedure Qualitative Test:

1. Bring reagents and specimens to room temperature before use.
2. Resuspend the Latex Reagent gently.
3. Place one drop (40 µl) of each test specimen or as reference one drop (40 µl) of the pos. / neg. control on to the reaction slide.
4. Add one drop Latex Reagent (40 µl) to each test field. Use spatula dispenser to spread reaction mixture over entire test field.
5. Rotate the slide for 2 minutes and read immediately under direct light.

Procedure Quantitative Test:

1. Bring reagents and specimens to room temperature before use.
2. Resuspend the Latex Reagent gently.
3. Using saline, dilute the specimens 1:2, 1 :4, 1 :8, 1 :16 or as needed.
4. Place one drop (40 µl) of each test specimen or as reference one drop (40 µl) of the pos. / neg. control on to the reaction slide.
5. Add one drop of the Latex Reagent (40 µl) to each test field. Use stick to spread reaction mixture over entire test field.
6. Rotate the slide for 2 minutes and read immediately under direct light.

Results Qualitative Test:

A negative reaction is indicated by a uniform milky suspension with no agglutination as observed with the negative control (2). A positive reaction is indicated by any observable agglutination in the reaction mixture. The specimen reaction should be compared with the positive control (1).

Results Quantitative Test:

The titer of the serum is the reciprocal of the highest dilution which exhibits a positive reaction.

Interpretation:

positive (1)

Agglutination
within 2 minutes

negative (2)

no Agglutination
within 2 minutes

Reference values:

Adults: < 8 IU/ml

The rheumatoid factors are immunoglobulins (most of the IgM) with antibody activity. These factors are present in most of patients suffering Rheumatoid Arthritis. There are different rheumatoid factors and does not exist any test capable to detect all of them, due to the fact that some of them act against human IgG, other against animal IgG, and other against both IgG. We recommend the use of WR test specific for detection of rheumatoid factors against animal IgG.

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range. For diagnostic purposes the rheumatoid factor results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

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