Optilene® Mesh LP

(en)	Instructions for use	– 2
(bg)	Инструкции за използване	– 4
(cs)	Návod k použití	– 6
(da)	Brugsvejledning	– 8
(de)	Gebrauchsanweisung	– 10
(el)	Οδηγίες χρήσης	– 12
(es)	Instrucciones de uso	– 14
(et)	Kasutusjuhend	– 16
(fi)	Käyttöohjeet	– 18
(fr)	Mode d'emploi	– 20
(hr)	Upute za uporabu	– 22
(hu)	Használati utasítás	– 24

(it)	Istruzioni per l'uso	– 26
(lt)	Naudojimo instrukcijos	– 28
(lv)	Lietošanas instrukcija	– 30
(nl)	Gebruiksaanwijzing	– 32
(no)	Bruksinstruksjoner	– 34
(pl)	Instrukcja stosowania	– 36
(pt)	Instruções de uso	– 38
(ro)	Instrucțiuni de utilizare	– 40
(sk)	Návod na použitie	– 42
(sl)	Navodila za uporabo	– 44
(sv)	Användarinstruktioner	– 46
(tr)	Kullanım talimatı	– 48

eIFU v.01 506677 - IFU(506371 01-07/25)

B | BRAUN

SHARING EXPERTISE

B. Braun Surgical, S.A.
Ctra de Terrassa, 121
08191 Rubí, Spain

Phone +34 93 586 62 00
Fax +34 93 588 10 96

www.bbraun.com

Optilene® Mesh LP

Description of the device

Optilene® Mesh LP is a synthetic non-bioabsorbable polypropylene surgical mesh which intended purpose is tissue reinforcement for preventing and repairing hernias.

Optilene® Mesh LP is a polypropylene mesh implant that has dyed blue stripes with copper phthalocyanine (Phthalocyaninato (2-) copper) to improve visibility and orientation control. It is characterized by its thin but extremely shape-stable mesh that after implantation adapts optimally to the longitudinal and latitudinal motion of the patient as a result of body movement and muscle contraction.

Indications for use

The indications for use of Optilene® Mesh LP are:

- Inguinal hernia repair
- Primary ventral hernia repair
- Incisional hernia repair
- Prevention of incisional hernia used as a prophylactic mesh in high-risk patients according to the professional criteria (such as but not limited to: age, nutrition, body mass index, comorbidities, etc.).

The target population for Optilene® Mesh LP is adult patients regardless of gender and who do not meet any contraindication.

Intended users are healthcare professionals like surgeons that are qualified and familiar with the surgical techniques.

Mode of action

Optilene® Mesh LP act as a reinforcement material in connective tissue weakness in order to repair or prevent hernia formation. In this sense, the mechanism of action is purely mechanical, acting as a barrier by distributing the internal pressure to the surface of the mesh, protecting the native fragile tissue from the pressure that causes or would potentially lead to hernia formation.

Contraindications

Optilene® Mesh LP should not be implanted:

- in contaminated and infected areas.
- in children, adolescents or patients who are still in the growth phase.
- in women with childbearing potential.
- in direct contact with the viscera (e.g.: intestines) in order to prevent adhesions.
- in pelvic floor.

Optilene® Mesh LP is contraindicated on patients with known severity or allergies to its components (refer to Materials Used, for details about components name)

Mode of application

Select Optilene® Mesh LP from the product range according to the size of the defect and cut it to fit the defect.

Optilene® Mesh LP is suitable for both conventional and laparoscopic application.

When repairing inguinal hernias, the mesh should, in order to prevent recurrence, be large enough to project beyond the pubic tubercle and pass reliably round the spermatic cord at the inner ring.

Meshes can be fixed with glues, sutures, staples, fibrin sealants or tack devices to be decided according to medical discretion. They should be placed at least 10 mm from the edge of the mesh.

Warnings

As with any implantable surgical device, strict aseptic techniques should be followed. If an infection develops, it should be treated aggressively. An unresolved infection may require removal of the device.

The presence of risk factors such as diabetes, hypertension or obesity might increase the probability of having follow-up complications.

Precautions

Before implanting Optilene® Mesh LP the surgeon must be familiar with the surgical technique, and specifically with the application of meshes and their fixation.

Opened unused or damaged packs must be discarded. Optilene® Mesh LP may not be used after the expiry date given on the packaging.

Do not reuse Optilene® Mesh LP, it is a single use product.

Do not reuse the product: infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

Of note, no special precautions are deemed necessary when performing MRI scans, since the implantable polypropylene meshes are non-metallic and non-conducting materials and therefore, considered Magnetic Resonance (MR) safe.

Side Effects

Side effects associated with the use of this device and common to other non-absorbable meshes may include local irritation in the area of the wound, inflammation due to foreign body reaction, seroma, hematoma, infection or fistula, delay in healing process, wound healing disorder, hemorrhages, ileus and intestinal lesion, migration, pain, chronic pain, feeling of foreign body, nerve damage, dysesthesia, organ prolapse and hernia recurrence. Systemic infections can lead to life threatening injuries.

Enterocutaneous fistula, bowel adhesions with consequent intestinal obstruction or perforation, chronic infection, sinus formation and/or erosion of the mesh into the viscera can occur as a result of direct contact of the mesh with the viscera.

Sterilization

Optilene® Mesh LP is sterilized by ethylene oxide and must not be resterilized. Only use Optilene® Mesh LP if the package is undamaged. The absence of channels along the sealing of the sterile barrier system must be checked as well as the absence of perforations on it. In case of identifying such defects, the device shall be discarded in the recommended manner.

Storage

No special storage conditions are required.

Device Disposal

Once the surgical procedure is completed, the different components of Optilene® Mesh LP have to be discarded in special intended containers. It is the responsibility of the user to determine if disposal material is hazardous according to federal, state and local regulations. Dispose of contents and container to comply with applicable local, state, national and international regulation.

Information to the User/Patient

The clinical staff shall complete the blank spaces of the supplied Implant Card (IC). The meaning of the symbols included in the IC is described in the Section "Symbol to be used on the label" of this Instructions for use. The IC shall be given to the patient after the surgical procedure.

In case that any serious incident occurs in relation to Optilene® Mesh LP, it should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Optilene® Mesh LP is permanent implantable mesh and therefore it is expected to last the entire life expectancy of the patient.

Optilene® Mesh LP is constructed with polypropylene undyed monofilament (~100 % w/w), and polypropylene dyed monofilament with colorant copper phthalocyanine (-2) (~100 % w/w of polypropylene and ≤ 0.5 % w/w of dyestuff) colorant in the blue stripes.

[<https://ec.europa.eu/tools/EUDAMED>]: Link to access at the Summary of Safety and Clinical Performance (SSCP) of Optilene® Mesh LP

Symbols to be used on the label

Do not re-use

Use-by date

Batch code

Do not re-sterilize

Polypropylene

Permanent undyed mesh with dyed stripes

Pore distribution (not true scaled)

Mesh geometry (not true scaled)

Sterilized using ethylene oxide

Single sterile barrier system with protective packaging inside

Do not use if package is damaged and consult instructions for use

Date of manufacture

eifu.bbraun.com

Users from the European Economic Area and Turkey must consult electronic Instructions for Use at eifu.bbraun.com. Users in other regions must consult the paper Instructions for Use supplied with the device.

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