Gazin®

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en Instructions for Use

Product description

Gazin® Dissecting Swabs, X-ray detectable, non-sterile, consist of 24-thread count cotton gauze in accordance with EN 14079. The X-ray detectable thread in the swabs consists of barium sulphate filaments and is wrapped with polypropylene.

Product composition

100% cotton; X-ray detectable thread (blue, polypropylene, polyester thread, barium sulphate)

Intended purpose

Gazin® gauze dissecting swabs are intended for the preparation and dissection of tissue layers and anatomical structures during surgical procedures.

Usage information/Precautions

This product is non-sterile. Gazin® Dissecting Swabs must be sterilised prior to surgically invasive use; resterilisation is not permitted. The product is intended for single use only.

Sterilisation instructions

1. Device requirements

Only sterilisers, which comply to the requirements of the norms EN 285 or EN 13060 and which use a pre-vacuum system for the sterilisation are allowed to be used. The use of the gravitational method is not allowed in the EU.

2. Sterilisation preparation

Before the product is sterilised, a suitable sterile barrier system must be selected, validated and used by the person responsible for the sterilisation (sterilising technician) in accordance with ISO 11607. The sterile barrier must be big enough for the product so that the seal is not under tension. The retail packaging is unsuitable for sterilisation.

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