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MARIDERM®

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Sterile, not resorbable suture material

Description

Mariderm® is a monofilament, not resorbable surgical suture material made of polyamide. Mariderm® fulfills the requirements of the European pharmacopeia and the United States pharmacopeia for sterile, not resorbable suture material, in the currently applicable versions. For better identification Mariderm® blue (Polyamide 6/66) is dyed with copper phthalocyanine blue (Colour-Index No. 71460) or Indigo Carmine (Colour-Index No. 73015) and Mariderm® black (Polyamide 6/6) with Hametin Logwood Black (Colour-Index No. 75290). Mariderm® is available in different thicknesses and lengths and as needle-suture-combination with atraumatic stainless steel needles of different types and sizes, as well as also without needle. Further details can be found in the catalogue.

Application areas

Mariderm® is indicated with applications where according to surgical practice not resorbable suture materials are used and where no permanent stability of the suture material is required. Mariderm® is largely suitable for all wound sutures and ligatures and is primarily indicated for skin suture. In addition to this Mariderm® black is indicated for microsurgery and ophthalmology. Apart from this polyamide sutures are often used for marking and as retaining sutures.

Contraindications

Mariderm® is not intended for application on the central cardio-vascular or nervous system. Due to tear strength decrease in the body Mariderm® must not be applied if a permanent stability of the suture material is required.

Mode of action

Not absorbable sutures cause minimal acute tissue reactions, which are followed by a slowly progressing encapsulation through connective tissue. Mariderm® is not resorbable, but can however due to hydrolysis of the polyamide in the body slowly decrease in tear strength. Therefore the application of Mariderm® is not suitable if a permanent stability of the suture material is required.

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Special warnings/ precautions for use/ interactions

With the application of absorbable suture material, the user should be familiar with surgical suturing methods and knotting techniques before using Mariderm®. The risk of wound dehiscence depends on the site of application and the applied suture material. Contaminated and infected wounds are to be provisioned according to generally valid surgical procedures. As with all foreign bodies a longer-term contact of the suture materials with saline solutions can lead to concrement formation (urinary or gall tracts).

During the handling damage to the suture threads is generally to be avoided. They must not be kinked or squeezed through surgical instruments such as pincette or needle holder. Needles should always be gripped between 1/3 and ½ of the length from the reinforced end to the needle tip. Deformation of the needle can lead to fracturing. The gripping of the needle at the needle tip can impair the penetration effect and lead to fracturing of the needle. Ionising radiation can have an effect on the properties of the material. With reuse hazard of infection for user and/or patient, as well as impairment of the functional ability of the suture material. Used needles must be professionally disposed of (for prevention of possible infection risks).

Adverse reactions

At the beginning low degree inflammatory tissue reactions can occur in the vicinity of the suture material. As with all foreign bodies not absorbable sutures can amplify a pre-existing infection.

Sterility

The suture material is sterile as long as the packing has not been opened or damaged. Damaged packs or packs appearing to be damaged are to be discarded. Unused suture material from opened packs must no longer be used! Do not re-sterilise!

Storage conditions

Do not store above 25°C and protect against moisture. Observe shelf-life!

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Stand 2019-12

Explanation of the symbols on the packaging

	Sterilised with ethylene oxide	
	Batch number	
	Article number	
	Observe instructions for use	
	Do not reuse	
	Use by	Year-Month (YYYY-MM)
	Manufacturer	Catgut GmbH Gewerbepark 18 08258 Markneukirchen Germany
	Do not use with damaged packing	
	Do not re-sterilise	
	Upper temperature limit value	
	Store dry	

	CE label with identification number of the designated authority	The product is consistent with the fundamental requirements of the directive 93/42/EEC related to medical devices.
--	---	--

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