BAUSCH + LOMB

CAPSULAR TENSION RING

ANILLO DE TENSIÓN CAPSULAR

KAPSELSPANNRING

ANNEAU DE TENSION CAPSULAIRE

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EN · ENGLISH

CAPSULAR TENSION RING

1 - IDENTIFICATION OF THE PRODUCT TO WHICH THE INSTRUCTIONS APPLY

Trade name	Reference	Description	UDI-DI Basic
Capsular Tension Ring	AC-PL	Round section ring preloaded	84351943010S2

2 - PRODUCT DESCRIPTION

The Capsular Tension Ring is a medical device that is implanted into the capsular bag to facilitate intraocular lens (IOL) implantation in patients with a weak or partially absent zonular body.

The Capsular Tension Ring is a one-piece, non-optical, round implant made entirely of white polymethylmethacrylate (PMMA). The composition of the material is as follows:

Compound	Concentration (% w/v)
Methyl methacrylate	99.69 - 99.71
Benzoyl peroxide	0.04 - 0.06
Tinuvin 326	0.25

The Capsular Tension Ring is a ring with an open area at a specific angle and two positioning holes.

3 - PRODUCT'S INTENDED PURPOSE

3.1. Intended use/specifications

The Capsular Tension Ring is designed to be inserted into the capsular bag at any time during cataract surgery to distend the capsular bag equatorially. The ring distributes traction forces on the capsule, expands the capsule, supports areas of zonular weakness, tightens the posterior capsule, redistributes tension in the capsular contracture syndrome and enters slightly off-centre capsular bags.

3.2. Indications

The Capsular Tension Ring is indicated for the stabilization of the capsular bag during and after cataract surgery in patients where it is present, among others:

- Primary or secondary zonular weakness
- Partial zonulolysis
- Iatrogenic Zonulopathy
- Eye trauma
- Pseudoexfoliation syndrome
- Risk of capsular retraction

3.3. Target patients

Capsular Tension Ring is indicated for adult patients undergoing cataract surgery where stabilisation of the capsular bag is required. There is no evidence on the safety of capsular tension rings during pregnancy and lactation in humans. Use during pregnancy and lactation is at the discretion of the ophthalmologist.

3.4. Intended users

The Capsular Tension Ring is designed for use only by ophthalmic surgeons who obtain the product directly from AJL or through an authorized distributor.

3.5. Lifetime

The lifetime of the Capsular Tension Ring is considered 15 years, as it is an implantable product intended for long-term use after implantation in the patient.

4 - CHOOSING THE RING

The diameter of the ring to be implanted must be chosen preoperatively, taking into account the White-white distance value of the eye. Equally, when implanting an intraocular lens, the geometry of the haptics must be taken into consideration when choosing the ring.

PRODUCT				Capsular Tension Ring		
GEOMETRY SECTION				Round section		
MODEL				AC-PL		
DIAMETER [mm]						
$\varnothing_o$	10,0	11,0	12,0	12,0	13,0	14,0
$\varnothing_c$	9,0	10,0	10,0	10,5	11,0	12,0

*PL= Preloaded; $\varnothing O$ = Open diameter; $\varnothing C$ = Closed diameter

5 - WARNINGS / CONTRAINDICATIONS

The surgeon must carry out a careful pre-operative evaluation and a clinical assessment to determine the risk-benefit ratio of

implantation in the following (non-exhaustive) list of pre-existing conditions:

- Zonular damage larger than 4 clock hours
- Microphthalmia
- Damage, tearing or perforation of the capsular bag

6 - SIDE EFFECTS AND RESIDUAL RISKS

As with any surgical procedure, cataract surgery with intraocular ring implantation, be it in conjunction with an intraocular lens or not, implies risks that the surgeon must evaluate as indicated in the relevant clinical literature. Postoperative complications may include, but are not limited to, the following:

- Dislocation of the capsular tension ring
- Dislocation of the intraocular lens
- Rupture of the Capsular Bag
- Glaucoma
- Iatrogenic zonular damage

The ring, injector or any part of it should not be re-sterilized, re-used and/or reprocessed as this may result in a number of known risks:

- Loss of product integrity and functionality.
- Potential changes in material and product design properties.
- Potential risk of infection or harm to the patient.

7 - PRECAUTIONS

1) Before surgery, the surgeon should consider the following:

- a. Do not store in direct sunlight or at temperatures 45°C.
- b. Do not use the injector if the level of pupillary dilation of the eye before surgery does not allow verification of the position of the capsular tension ring within the capsular bag after surgery.

2) Do not use the product in case of:

- a. The expiration date has expired.
- b. Lack of labelling and/or instructions for use.
- c. Unintentional opening of the package before the intended date of use. The sterility of the contents can only be guaranteed if the protective packaging of the product has not been opened or damaged.
- d. If the injector has been dropped or if any part of the injector has been accidentally knocked or damaged while out of its protective packaging.


3) During surgery, the surgeon should consider the following

- a. Do not touch the injector tip to avoid damage.
- b. If the injector requires excessive force to eject the ring, discard the product.
- c. Do not implant in the sulcus or anterior chamber, as the intended implantation site is the capsular bag.

8 - INSTRUCTIONS FOR USE / PRODUCT HANDLING

Check that the expiry date of the product has not expired. Once you have opened the box, check that the label on the container matches the model and diameter that you need and that it is not damaged.

1. In aseptic conditions, open the blister slowly and continuously and carefully remove the injector from the blister (Figure 1).
2. Turn the wheel forward until the ring is fully inserted into the injector (Figure 2).
3. Unlock the injector by sliding the tab towards the centre until you hear a “click” (Figure 3).
4. Insert the tip of the injector into the capsular bag (Figure 4).
5. Turn the wheel backwards to release the ring in a smooth and controlled manner along the capsular equator (Figure 5).
6. The ring will free itself from the hook when it is fully released from the injector.
7. Turn the wheel forward until the hook is retracted and remove the injector from the eye.

5. Rotate the injector to adjust the output of the AC/ACS rings

9 - PRODUCT DISPOSAL

It must be disposed of in accordance with local regulatory requirements. Discarded rings (whether used or unused) are classified as waste.

10 - CLINICAL BENEFITS

Implantation of the Capsular Tension Ring may be beneficial in the following (non-exhaustive) list of situations:

- Stabilize and centre the capsular bag
- To promote centration and stabilization of intraocular lenses

11 - INSTRUCTIONS FOR COMPLETING THE IMPLANT CARD

Capsular Tension Ring is accompanied by an implant card that must be completed by the healthcare professional and delivered to the patient as a permanent record of their implant; so that they can present it to any ophthalmologist, they may consult in the future. It must include information on the patient, the doctor and the health centre where the surgical procedure was carried out, as well as data on the implant itself by means of the label provided for that purpose.

Below are instructions on how to complete the implant card:

- 1) Attach the label provided in the secondary packaging. This label is identified by the phrase "Stick on the implant card" at the top.
- 2) Select the eye in which the Capsular Tension Ring has been implanted.
- 3) Fill in the patient's name or ID.
- 4) Write down the date of implantation.
- 5) Indicate the address of the healthcare professional or healthcare institution where the patient information can be found.


Figure 6. Back of implant card


Figure 7. Front of implant card

The healthcare professional must inform the patient of the importance of carrying the implant card with them and the possibility of accessing the patient information sheet from the website mentioned on the implant card (www.patientinfo.ajlsa.com), as well as their right to request such information in physical format if they do not have internet access. In the latter case, the doctor must provide the patient with a physical copy

of the patient information sheet, which they can access via the above-mentioned link. If this is not possible, they must contact the manufacturer to obtain a physical copy.

12 - NOTIFICATION

Any adverse events or complaints should be reported to AJL's Quality Control Department by email to ajlsa@ajlsa.com. In the event of a serious incident, the competent national authority should be notified.

13 - CLINICAL SAFETY AND PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) of products registered under Regulation (EU) 2017/745 will be published in the European Medical Device Database (EUDAMED) once it is operational.

The address for accessing the product's SSCP is <https://ec.europa.eu/tools/eudamed/#/screen/home> and it can be accessed with the basic UDI-DI code corresponding to the product, which is 84351943010S2.

14 - LIMITATION OF WARRANTY AND LIABILITY

Placement of the Capsular Tension Ring is a surgical procedure and is associated with some risks. AJL Ophthalmic S.A. has provided the necessary information regarding these risks and on the methods and techniques to be used. The doctor must provide the patient with all relevant information. In particular, AJL Ophthalmic S.A. is exempt from any liability related to any harm and injury to the patient associated with:

- The implantation method or procedure used by the doctor should they fail to follow the recommendations provided by AJL Ophthalmic S.A.
- Incorrect prescription, selection or use of the Capsular Tension Ring for a particular patient.

For exchanges or returns, please contact the manufacturer or local distributor. In the event that the product is damaged or non-sterile, it must be returned in its original packaging bearing the serial number, reference, and reason for the return, to AJL Ophthalmic S.A. or to the local distributor.

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Medilas AG,
AJL OPHTHALMIC, S.A.