RÜSCH Tracheal tubes made of PVC for single use Instructions for use	en
RÜSCH Einmal-Trachealtuben aus PVC Gebrauchsanleitung	de
RÜSCH Tubes trachéaux á usage unique en PVC Mode d'emploi	fr
RÜSCH Tubos endotraqueales desechables de PVC Instrucciones de uso	es
RÜSCH Tubi tracheali monouso in PVC Istruzioni d'uso	it
RÜSCH Tubos endotraqueais de utilização única feitos de PVC Instruções de uso	pt
RÜSCH PVC-trakealtuber för engångsbruk Bruksanvisning	sv
RÜSCH Tracheale tubes van PVC voor eenmalig gebruik Gebruiksaanwijzingen	nl
RÜSCH Tracheal tuber til éngangs-brug Brugsanvisning	da
RÜSCH-trakeaaliputket valmistettu PVC:stä, kertakäyttöön Käyttöohje	fi
RÜSCH Luftrørsslanger laget av PVC, til engangsbruk Bruksanvisning	no
RÜSCH Tek kullanımlık PVC trakea tüpleri Kullanım bilgileri	tr
RÜSCH Jednorazowych rurek dotchawiczych firmy wykonanych z PVC Instrukcja użycia	pl
RÜSCH Τραχειακοί σωλήνες από PVC μίας χρήσης Οδηγίες χρήσης	el
RÜSCH Эндотрахеальная трубка с манжетой и портом для аспирации надманжеточного пространства Инструкции по применению	ru
Trahealni tubusi RÜSCH izdelani iz PVC za enkratno uporabo Navodila za uporabo	sl
Tracheální kanyly RÜSCH vyrobené z PVC, na jedno použití Návod k použití	cs
Tracheálne trubice RÜSCH vyrobené z materiálu PVC na jednorazové použitie Návod na použitie	sk
Трахеални тръби RÜSCH, изработени от PVC за еднократна употреба Инструкции за употреба	bg
Trahealne cijevi RÜSCH od PVC-a, za jednokratnu uporabu Upute za uporabu	hr
Intubatsioonitorud RÜSCH, valmistatud PVC-st, ühekordseks kasutamiseks Kasutusjuhend	et
RÜSCH légcső kanülök PVC-ből készült, egyszeri használatra Használati utasítás	hu
Tuburi traheale RÜSCH fabricate din PVC, de unică folosință Instrucțiuni de utilizare	ro

3 Instructions for use	en	☐
5 Gebrauchsanleitung	de	☐
7 Mode d'emploi	fr	☐
9 Instrucciones de uso	es	☐
11 Istruzioni d'uso	it	☐
13 Instruções de uso	pt	☐
15 Bruksanvisning	sv	☐
17 Gebruiksaanwijzingen	nl	☐
19 Brugsanvisning	da	☐
21 Käyttöohje	fi	☐
23 Bruksanvisning	no	☐
25 Kullanım bilgileri	tr	☐
27 Instrukcja użycia	pl	☐
29 Οδηγίες χρήσης	el	☐
31 Инструкции по применению	ru	☐
33 Navodila za uporabo	sl	☐
35 Návod k použití	cs	☐
37 Návod na použitie	sk	☐
39 Инструкции за употреба	bg	☐
41 Upute za uporabu	hr	☐
43 Kasutusjuhend	et	☐
45 Használati utasítás	hu	☐
47 Instrucțiuni de utilizare	ro	☐

Important!

Read the entire instructions prior to use.

This instruction insert applies to the following Rüsch products:

Standard tracheal tubes, pre-formed tracheal tubes, armoured tracheal tubes, special tracheal tubes, tracheal tubes with and without cuff, tracheal tubes with and without intubation stylet.

Additional Rüsch products, which are added after printing of these instructions for use, will be labelled separately.

Description:

Rüsch tracheal tubes made of polyvinyl chloride (PVC) can be supplied depending on the model as standard tracheal tubes, tracheal tubes with a metal spiral embedded in the tube wall, and with or without cuffs respectively. Depending on the model, the filling system consists of a filling tube, a pilot balloon, a self-closing valve with a Luer connection or a universal syringe tip.

Depending on the model, Rüsch tracheal tubes are available

- with or without Murphy eye
- as a nasal/oral
- oral or nasal version
- and in a curved, anatomically shaped or straight form.

Not all versions are available in all models and sizes.

Note:

All of the information relating to the cuffed versions does not apply to tracheal tubes without cuff.

Indications:

Rüsch tracheal tubes are indicated for oral or nasal intubation for airway management. Reinforced tracheal tubes may be used to reduce the potential for kinking wherever unusual positioning of the head and neck is required following intubation.

Contraindications:

None known.

Adverse reactions:

Reported adverse reactions associated with tracheal tubes are many and diverse.

Standard textbooks should be consulted for specific adverse reaction information (e.g. J.L Benumof, Airway Management, Mosby, 1996).

SUGGESTED DIRECTIONS FOR USE

1. Carefully remove the sterile tracheal tube from its protective packaging.
2. Firmly seat the 15 mm connector into the tracheal tube, to help prevent disconnection during use. (Not necessary where the connector is already firmly seated or glued into the tube)
3. The patency of the main tube lumen should be checked prior to intubation.
4. If the tube has to be cut to length prior to intubation, remove the fitted 15 mm connector and cut the tube at the appropriate point. The 15 mm connector should then be firmly reinserted into the tracheal tube.
Consideration should be given to the specific WARNINGS/PRECAUTIONS (general) in the instruction insert.
5. Check the cuff, pilot balloon and valve or other filling systems on its leak-proofness and integrity prior to use. Insert a Luer syringe into the Luer adapter for the filling system and inject enough air to fully inflate the cuff.
6. After the test inflation of the cuff, completely evacuate the air.
7. Intubate the patient following currently accepted medical techniques noting the specific WARNINGS/PRECAUTIONS (cuff-related) in the instruction insert.
8. Once the patient is intubated, use a Luer syringe or an appropriate cuff pressure gauge (e.g. Rüsch ENDOTEST) to inflate the cuff with just enough air to provide an effective seal of the trachea. The use of minimal occluding volume, minimum leak techniques and continuous monitoring of the cuff pressure can help reduce the incidence of many of the adverse reactions associated with the use of cuffed tracheal tubes.
9. Fix the tube to the patient following currently accepted medical techniques.
10. In the case of tracheal tubes that include a valve, remove the syringe from the valve after cuff inflation. Leaving the syringe attached will keep the valve open, permitting the cuff to deflate.
11. Check that the inflation system is not leaking. Cuff pressure should be measured periodically or permanently. The integrity of the system should be verified periodically during the intubation period. Any deviation from the selected sealing pressure should be investigated and corrected immediately.
12. Prior to extubation, deflate the cuff by inserting a syringe into the cuff filling system and removing all of the air contained within the cuff until the pilot balloon has collapsed.
13. Extubate the patient following currently accepted medical standards.
Dispose of the tracheal tube. The specific WARNING/PRECAUTIONS (general) in the instruction insert should be noted.

WARNINGS / PRECAUTIONS (cuff-related)

- Inflation of the cuff by "feel" alone or by using a measured amount of air is not recommended, since resistance is an unreliable guide during inflation. Diffusion of nitrous oxide mixture, oxygen or air may either increase or decrease the cuff volume and the cuff pressure. The use of a cuff pressure gauge (Rüsch ENDOTEST, REF 112700) will help to monitor and adjust the cuff pressure.
- Do not overinflate the cuff. The cuff pressure should normally not exceed 25 cm H₂O. Overinflation can result in tracheal damage, rupture of the cuff with subsequent deflation, or cuff distortion, which may lead to airway blockage.
- The use of lidocaine topical aerosols has been associated with the formation of pinholes in PVC cuffs (Jayasuiya KD, Watson WF : "P.V.C. cuffs and lidocaine-based aerosol"; Br J Anaesth. 1981 Dec; 53 (12): 1368). Expert clinical judgement is required when using the substance, to help prevent cuff leaks. The use of a water soluble lubricant gel is recommended.

- Various bony anatomical structures (e.g. teeth, turbinates) within the intubation routes or any intubation aid with sharp surfaces can damage the integrity of the cuff. Care must be taken to avoid damaging the thin-walled cuff during intubation, which could result in the patient requiring the trauma of extubation and re-intubation. If the cuff is damaged, the tube should not be used.
- Deflate the cuff prior to repositioning the tube. Movement of the tube with the cuff inflated could result in damage to the cuff or in patient injury, requiring possible medical interventions.
- Syringes, stopcocks or other devices should not be left inserted in the inflation system for extended periods of time.

WARNINGS / PRECAUTIONS (general)

- Using tracheal tubes during procedures which will involve the use of a laser or an electrocautery with the tracheal tube, especially in the presence of oxygen-enriched or nitrous oxide-containing mixtures, could result in rapid combustion of the tube with harmful thermal effects and with emission of corrosive and toxic combustion products including hydrochloric acid (HCL).
- In case the tube is being cut to length, when reinserting the 15 mm connector, ensure that it is not reinserted in the proximity of the joint between the inflation system and the tube. If this occurs, accidental blockage or damage with subsequent leakage of the inflation system could result.
- RÜSCH armoured tracheal tubes or tubes with a firmly glued 15 mm connector cannot be cut to length.
- The user should be aware of anatomical variations including the length of the airway. Reliance on the cm graduation markers, black points or tips on tracheal tubes or black circular marks above the cuff should never be a substitute for expert clinical judgement.
- When the position of the tube or of the patient himself is altered after intubation, it is essential to verify that the tube remains correctly positioned.
- If intubation is performed using a stylet, verify that the end of the stylet does not protrude from the end of the tube or pass through the Murphy eye (if present).
- Should extreme flexion (chin-to-chest) of the head or movement of the patient (e.g. to the lateral or prone position) be anticipated after intubation, use of RÜSCH armoured tracheal tubes should be considered.
- The use of an appropriate bite block is still required even when armoured tracheal tubes are used.
- Expert clinical judgement should be exercised in the selection of the appropriate size of tracheal tube for each individual patient.
- Intubation and extubation should be performed following currently accepted medical techniques.
- If the tracheal tube is lubricated prior to intubation, it is essential to verify that the lubricant does not enter and occlude the tube lumen or the cuff filling system, thereby preventing ventilation, or causing damage to the cuff.
- Use of lubricating solutions to ease reinsertion of the 15mm-connector is not recommended as it may contribute to accidental disconnection.
- **WARNING :** Do not use lubricants with a petroleum base.
- The device should be disposed of in accordance with applicable national regulations.
-  Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
- Some tracheal tubes contain phthalates : See product label

Contains DEHP

The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction.

Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application.

Medical devices containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Delivery:

Sterile unless packaging has been opened or damaged.

Discard after single use, do not re-use

Sterile, see packaging.

Exposure to elevated temperature and UV light and excessive moisture should be avoided during storage.

Storage Instruction:

Keep away from sunlight and keep dry. Do not use if the product sterilisation barrier or its packaging is compromised.

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