5

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Front thigh muscles

Sitting. The back is straight, both legs are set up, toe tips are pulled up. SISSEL® Fit-Tube centered around both feet, hands grab loops. Both hands rest on the ground next to the body at pelvic height. Stretch legs alternately against the tube resistance just above the ground. The tip of the toe remains pulled up.

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Buttock muscle

Standing position. Shift weight to the front leg, slightly bent leg. Front foot stands on SISSEL® Fit-Tube, hand of the same side grabs loop. Upper arm and elbow are attached to the body. Place the second loop of the SISSEL® Fit-Tube around the foot or ankle of the rear leg. Upper body upright, lower shoulders to the rear-bottom. Back leg stretched slowly to the back and forward again. Keep your foot stretched

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Upper arm: Biceps

Standing Position. Shift weight to the front leg, both knees slightly bent. Front foot stands on SISSEL® Fit-Tube, hands grab loops. Upper arms and elbows are attached to the body. Upper body upright, lower shoulders to the rear-bottom. Alternately slow bending and stretching the elbows.

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Upper arm: Triceps

Standing position. Shift weight to the front leg, both knees slightly bent. Front foot stands on SISSEL® Fit-Tube, hands grab loops. Arms long, slightly bent in the elbow, next to the body. Upper body upright, lower shoulders to the rear-bottom. Alternately lead the arms backwards and forwards again.

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Straight abdominal muscles

Standing position. Both legs are set up, toe tips are pulled up. SISSEL® Fit-Tube lies slightly stretched under the shoulder blades. Hands grab loops, arms stretched up. Lift the head and upper body slightly under exhalation, the stretched arms lead straight towards the knees. Go back with the inhalation.

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Lateral abdominal muscles

Supine position. Both legs are set up, toe tips are pulled up. SISSEL® Fit-Tube lies slightly stretched under the shoulder blades. Hands grab loops, arms stretched up. Lift the head and upper body slightly under exhalation, the stretched arms lead sideways towards the knees. Go back with the inhalation. Left and right alternately.

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Precautionary instructions

- Check integrity of the product, do not use any damaged product.
- Do not apply product in the immediate vicinity of heat sources, protect from direct sunlight.
- Keep product away from sharp objects or edges.
- Use during pregnancy as well as health problems after medical consultation.
- The product shall not be kept outside the reach of children. Danger of strife.

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SISSEL® GmbH · Bruchstrasse 48 · 67098 Bad Dürkheim · Germany