Optilene® Mesh

(en)	Instructions for use	– 2
(bg)	Инструкции за използване	– 4
(cs)	Návod k použití	– 6
(da)	Brugsvejledning	– 8
(de)	Gebrauchsanweisung	– 10
(el)	Οδηγίες χρήσης	– 12
(es)	Instrucciones de uso	– 14
(et)	Kasutusjuhend	– 16
(fi)	Käyttöohjeet	– 18
(fr)	Mode d'emploi	– 20
(hr)	Upute za uporabu	– 22
(hu)	Használati utasítás	– 24

(it)	Istruzioni per l'uso	– 26
(lt)	Naudojimo instrukcijos	– 28
(lv)	Lietošanas instrukcija	– 30
(nl)	Gebruiksaanwijzing	– 32
(no)	Bruksinstruksjoner	– 34
(pl)	Instrukcja stosowania	– 36
(pt)	Instruções de uso	– 38
(ro)	Instrucțiuni de utilizare	– 40
(sk)	Návod na použitie	– 42
(sl)	Navodila za uporabo	– 44
(sv)	Användarinstruktioner	– 46
(tr)	Kullanım talimatı	– 48

CE 0123

eIFU v.01 506679 - IFU(506368 01-06/25)

B | BRAUN

SHARING EXPERTISE

B. Braun Surgical, S.A.
Ctra de Terrassa, 121
08191 Rubí, Spain

Phone +34 93 586 62 00
Fax +34 93 588 10 96

www.bbraun.com

Optilene® Mesh

Description of the device

Optilene® Mesh is a synthetic non-bioabsorbable polypropylene surgical mesh which intended purpose is tissue reinforcement for repairing hernias.

Optilene® Mesh is a polypropylene mesh implant that has dyed blue stripes with copper phthalocyanine (Phthalocyaninato (2-) copper) to improve visibility and orientation control. It is characterized by its thin but extremely shape-stable mesh that after implantation adapts optimally to the longitudinal and latitudinal motion of the patient as a result of body movement and muscle contraction.

Indications for use

The indication for use of Optilene® Mesh are:

- Inguinal hernia repair
- Primary ventral hernia repair
- Incisional hernia repair

The target population for Optilene® Mesh is adult patients regardless of gender and who do not meet any contraindication.

Intended users are healthcare professionals like surgeons that are qualified and familiar with the surgical techniques.

Mode of action

Optilene® Mesh act as a reinforcement material in connective tissue weakness in order to repair hernia. In this sense, the mechanism of action is purely mechanical, acting as a barrier by distributing the internal pressure to the surface of the mesh, protecting the native fragile tissue from the pressure that causes or would potentially lead to hernia formation.

Contraindications

Optilene® Mesh should not be implanted:

- in contaminated and infected areas.
 - in children, adolescents or patients who are still in the growth phase.
 - in women with childbearing potential.
 - in direct contact with the viscera (e.g.: intestines) in order to prevent adhesions.
 - in pelvic floor.
- Optilene® Mesh are contraindicated on patients with known severity or allergies to its components (refer to Materials Used, for details about components name).

Mode of application

Select Optilene® Mesh from the product range according to the size of the defect and cut it to fit the defect.

Optilene® Mesh is suitable for both conventional and laparoscopic application. When repairing inguinal hernias, the mesh should, in order to prevent recurrence, be large enough to project beyond the pubic tubercle and pass reliably round the spermatic cord at the inner ring.

Meshes can be fixed with glues, sutures, staples, fibrin sealants or tack devices to be decided according to medical discretion. They should be placed at least 10 mm from the edge of the mesh.

Warnings

As with any implantable surgical device, strict aseptic techniques should be followed. If an infection develops, it should be treated aggressively. An unresolved infection may require removal of the device.

The presence of risk factors such as diabetes, hypertension or obesity might increase the probability of having follow-up complications.

Precautions

Before implanting Optilene® Mesh the surgeon must be familiar with the surgical technique, and specifically with the application of meshes and their fixation.

Opened unused or damaged packs must be discarded. Optilene® Mesh may not be used after the expiry date given on the packaging.

Do not reuse Optilene® Mesh, it is single use product.

Do not reuse the product: Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

Of note, no special precautions are deemed necessary when performing MRI scans, since the implantable polypropylene meshes are non-metallic and non-conducting materials and therefore, considered Magnetic Resonance (MR) safe.

Side Effects

Side effects associated with the use of this device and common to other non-absorbable meshes may include local irritation in the area of the wound, inflammation due to foreign body reaction, seroma, hematoma, infection or fistula, delay in healing process, wound healing disorder, hemorrhages, ileus and intestinal lesion, migration, pain, chronic pain, feeling of foreign body, nerve damage, dysesthesia, organ prolapsed and hernia recurrence. Systemic infections can lead to life threatening injuries.

Enterocutaneous fistula, bowel adhesions with consequent intestinal obstruction or perforation, chronic infection, sinus formation and/or erosion of the mesh into the viscera can occur as a result of direct contact of the mesh with the viscera.

Sterilization

Optilene® Mesh is sterilized by ethylene oxide and must not be resterilized.

Only use Optilene® Mesh if the package is undamaged. The absence of channels along the sealing of the sterile barrier system must be checked as well as the absence of perforations on it. In case of identifying such defects, the device shall be discarded in the recommended manner.

Storage

No special storage conditions are required.

Device Disposal

Once the surgical procedure is completed, the different components of Optilene® Mesh have to be discarded in special intended containers. It is the responsibility of the user to determine if disposal material is hazardous according to federal, state and local regulations. Dispose of contents and container to comply with applicable local, state, national and international regulation.

Information to the User/Patient

The clinical staff shall complete the blank spaces of the supplied Implant Card (IC). The meaning of the symbols included in the IC is described in the Section "Symbol to be used on the label" of this Instructions for use. The IC shall be given to the patient after the surgical procedure.

In case that any serious incident occurs in relation to Optilene® Mesh, it should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Optilene® Mesh are permanent implantable meshes and therefore they are expected to last the entire life expectancy of the patient.

Optilene® Mesh are constructed with polypropylene undyed monofilament (~100 % w/w), and polypropylene dyed monofilament with colorant copper Phthalocyaninate (-2) (~100 % w/w of polypropylene and ≤ 0.5 % w/w of dyestuff) in the blue stripes.

[<https://ec.europa.eu/tools/EUDAMED>]: Link to access at the Summary of Safety and Clinical Performance (SSCP) of Optilene Mesh

Symbols to be used on the label

	Do not re-use
	Use-by date
	Batch code
	Do not re-sterilize
	Polypropylene
	Permanent undyed mesh with dyed stripes
	Pore distribution (not true scaled)
	
	Mesh geometry (not true scaled)
	Sterilized using ethylene oxide
	Single sterile barrier system with protective packaging inside
	Do not use if package is damaged and consult instructions for use
	Date of manufacture
	Users from the European Economic Area and Turkey must consult electronic Instructions for Use at eifu.bbraun.com . Users in other regions must consult the paper Instructions for Use supplied with the device.
	Medical device
	Catalogue number
	Unique Device Identifier
	Patient identification
	Date
	Health care centre or doctor
	Patient information website
	Manufacturer

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