erbe

30195-011_V24946
2024-10
Sheet 1 / 6

Manufacturer

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Bipolar forceps

1 General instructions for use

These notes on use describe the normal use of the product.

Please read all information carefully.

These notes on use do not replace the user manuals of the units and the products that are used with the instrument! Read the user manuals of the units and products and ask Erbe or your distributor in case of doubt!

Note: Report serious incidents with the product to your local dealer or Erbe. If you are a user in the European Union, also report incidents to the responsible authority in your Member State.

1.1 Intended use / Indications for use

The bipolar forceps are intended for coagulation of tissue in open surgical procedures.

1.2 Clinical benefit

- hemostatic coagulation

1.3 Compatibility

When attaching the bipolar forceps to electrosurgical units, use the following Erbe connecting cables for bipolar forceps: 20196-045, 20196-047, 20196-053, 20196-055, 20196-057, 20196-059, 20196-061.

4 How to use

Clean, disinfect, and sterilize the product before using it for the first time and before every subsequent use.

4.1 Checking the product

WARNING! Check the product for external damage before every use. You must observe the notes in the chapter “Check”. Do not use this product if it is damaged.

4.2 Connecting the product

WARNING! Risk of injury! Do not activate the electrosurgical unit while you are connecting the instrument.

- Attach the product to a matching Erbe connecting cable for bipolar instruments.
- Connect the product to the bipolar socket of the electrosurgical unit.

4.3 Selecting unit settings

WARNING! Risk of accidental tissue damage. Always use the lowest power and voltage settings for the required surgical effect.

Only select unit settings that do not exceed the maximum electrical capacity of the product (V_p max or kV_p max) as the product may otherwise be damaged. Notes on selecting appropriate unit settings can be found in the user manual of the electrosurgical unit (for Erbe electrosurgical units, read the "Accessories" section).

Erbe electrosurgical unit: Erbe recommends operating the instrument in the soft coagulation mode. Sparks generated in the forced coagulation mode may cause increased wear.

4.4 Conducting a performance test

Carry out the performance test described here before every application.

- Moisten a gauze strip with a sterile NaCl solution.
- Grasp the gauze strip with the instrument jaws and activate the HF current.

Steam will form after approx. 2 to 3seconds if the instrument is working properly.

Do not use this product if it is damaged!

4.5 During use

CAUTION! Deposits may impair the effectiveness of the instrument. Regularly clean the instrument tips of any blood and tissue residues.

4.6 Precleaning the product in the procedure room

Use a soft cloth (such as a compress) to remove any obvious contamination immediately after use and/or place the product in deionized water.

4.7 After use

Start processing the product within 20 minutes of use.

5 Cleaning, disinfection, sterilization

5.1 Safety instructions for processing

Suitable for cleaning in an ultrasonic bath (max. 10 minutes per processing cycle).

On no account use sharp or abrasive objects for cleaning purposes.

Disinfectants must be rinsed off thoroughly after use.

Maximum air pressure for drying: 2 bar.

Maximum water pressure for rinsing: 2 bar.

During machine cleaning and disinfection, do not exceed 93 °C.

Observe the specifications of the manufacturer of the detergent and disinfectant and the information in the chapter “Overview of validated procedures”.

During sterilization, do not exceed 134 °C.

Observe the national regulations for processing products after treating patients with Creutzfeldt-Jakob disease (CJD), with suspected CJD, or potential variants of CJD.

Observe the national regulations for processing parameters.

Do not sterilize in hot air.

Do not sterilize with gamma radiation.

Do not sterilize with STERRAD (J&J)!

5.2 Processing limitation

Frequent processing and surgical stress will affect this product. Processing and type of application may generally affect the service life of the product. Follow the safety instructions when inspecting the product. Do not continue to use the product in case of obvious damage or malfunction.

Detailed visual inspection and function testing before the next use is the best way to identify a product that is no longer functional; see Check.

5.3 Dismantling

- Disconnect the connecting cable from the product.

5.4 Required aids

Processing step	Items needed
Precleaning	Soft plastic brush/soft cloth
Manual cleaning/ disinfection	Soft plastic brush/soft disposable cloth (low particle) Compressed air gun (2 bar maximum)

5.5 Precleaning

WARNING! Risk of protein fixation. Avoid protein-fixing conditions, e.g., high temperatures. Use cold water (< 55 °C).

Use water of at least drinking quality for precleaning and, if necessary, a non-fixing disinfectant.

Suitable for cleaning in a 3% hydrogen peroxide bath (max. 30 minutes per processing cycle).

- Remove surface contaminants with a soft brush/cloth while holding the product in a water bath and/or rinsing the product under running water.

5.6 Manual cleaning and disinfection

Use water of at least drinking quality, and for the final rinse use sterile deionized water.

Use a liquid detergent suitable for the preparation of a dipping bath.

Use a disinfectant compatible with the detergent suitable for the preparation of a dipping bath.

The detergent and disinfectant must be suitable for medical devices made of plastic and metal and have a pH value between 5.5 and 12.3.

Observe the specifications from the manufacturer of the detergent and disinfectant and in the chapter “Overview of validated procedures”.

- Prepare a cleaning bath and a separate disinfectant bath.
- Immerse the product completely in a cleaning bath, ensuring that the product does not come into contact with other components in the bath.
- Thoroughly clean the surfaces with a soft brush/soft disposable cloth.
- Rinse the surfaces with sterile deionized water (at least 1 minute).
- Repeat the preceding cleaning steps if the product still shows visible contamination.
- Immerse the product completely in a disinfectant bath, ensuring that the product does not touch other components in the bath.
- Rinse the surfaces with sterile deionized water (at least 1 minute).
- Dry the product until there are no more fluid residues visible (e.g. with filtered compressed air).

5.7 Cleaning and disinfection by machine

Use water of at least drinking quality, and for the final rinse use deionized water (microbiologically at least drinking water quality).

The washer/disinfector must have fundamentally certified effectiveness (e.g. CE marking according to DINENISO15883).

The detergent and disinfectant must be suitable for medical devices made of plastic and metal and have a pH value between 5.5 and 12.3.

Observe the specifications from the manufacturer of the detergent and disinfectant and in the chapter “Overview of validated procedures”.

- Place the product in a suitable rinsing basket carefully. In doing so make sure the product is not touching other instruments or parts of instruments.
- Start a tested program with the following properties:
 - Sufficient cleaning effect.
 - Thermal disinfection: 5 to 10 minutes at 90 to 93 °C, $A_0 \geq 3000$.
 - Adequate rinse cycles
 - Final rinse with deionized water (microbiologically at least drinking water quality).
 - Sufficient product drying (no visible fluid residues).
- If contaminants are still visible on the product after the end of the routine, repeat precleaning and machine cleaning and disinfection.

5.8 Check

WARNING! Check the product for visible damage and wear. You can perform the visual inspection with or without magnification.

- Check the following:
 - Damage to the product, e.g. cracks, rough surface, splitting, corrosion.
 - Damage to the product's insulation and/or the cable/plug, e.g. cracks and fractures. Since under certain circumstances it may not be possible to detect microscopic cracks in the insulation, even under magnification, Erbe recommends also testing the insulation with an electric testing device.
 - Bending of the product or parts of the product. Do not attempt to straighten bent instruments!
 - Legibility of the inscription

Discoloration can be removed by thoroughly cleaning with a cleaning agent containing phosphoric acid. Only the instrument tip may come into contact with the cleaning solution.

Do not use this product if it is damaged!

5.9 Packaging

- Pack the product in disposable sterilization packaging (single or double packaging) made from paper/foil and/or in a sterilization container as specified in the EN ISO 11607 series.

5.10 Sterilization

Only sterilize products that have been cleaned and disinfected.

Erbe Elektromedizin recommends steam sterilization with the method described below. If other sterilization methods are used, Erbe Elektromedizin shall bear no responsibility.

Steam sterilization

- Fractionated vacuum method with adequate product drying
- Holding time 3 to 20 minutes at 132 to 134 °C
- Sterilizers in accordance with applicable national standards and regulations (e.g., DINEN13060 or DINEN285)
- Sterilization process validated according to DINEN ISO17665

Please follow the recommendations issued by the sterilizer manufacturer concerning loading, handling, and drying times.

5.11 Overview of validated procedures

The following procedures have been validated as suitable for processing the product:

- Manual cleaning/disinfection with the detergent Cidezyme and disinfectant Cidex OPA (Johnson & Johnson Medical Limited, Gargrave, Skipton).
 - Cidezyme: Dosage 1.6%, 5 minutes exposure time at room temperature.
 - Cidex OPA: Dosage 0.3%, 12 minutes exposure time at 20 °C.
- Machine cleaning and disinfection in a G 7836 CD disinfector (Miele & Cie. KG, Gütersloh, Germany) with the cleaning agent neodisher mediclean forte (Dr. Weigert GmbH & Co. KG, Hamburg, Germany).
 - neodisher mediclean forte: Dosage 0.5%, 5 minutes exposure time at 55 °C.
 - Rinsing: 2 minutes at 10 °C
 - Disinfection: 5 minutes at 90 °C.
- Sterilization in an autoclave Systec V-150 (SystecGmbH, Linden/Germany).
 - Steam sterilization with saturated steam, fractionated vacuum procedure, 3 minutes, 132 °C, 3 bar.

Erbe recommends the processing procedures described in these notes on use. Equivalent different procedures are possible if not explicitly excluded. It is incumbent on the user to ensure the suitability of the actual procedures used by means of suitable measures (e.g., validation, routine monitoring, check of material compatibility).

6 Disposal

WARNING! Risk of infection by contaminated and potentially infectious material.

Dispose of the product, packaging material, and accessories (if available) in accordance with the national guidelines and legislation applicable in each case.

7 Symbols

Symbol	Explanation
	Consult instructions for use
	Caution, read safety instructions before use.
	Catalogue number
	Batch code
	Medical device
	Manufacturer
	Date of manufacture
	Keep away from sunlight
	Keep dry
	Quantity (x)
	European conformity marking

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© Erbe Elektromedizin GmbH
- Systec V-150 otoklavda sterilizasyon (Systec GmbH, Linden/Almanya).