Tegaderm™ Alginate			GB
High Integrity Alginate Dressings			
Description: Tegaderm™ Alginate dressing is a sterile, non-woven pad or tape composed of calcium/sodium alginate. Tegaderm™ Alginate dressing is soft, highly absorbent, conformable and fast gelling. Not made with natural rubber latex.			
Indications: Tegaderm™ Alginate dressing is indicated for the following medical conditions to heavily exuding chronic and acute wounds: pressure ulcers, venous and arterial leg ulcers, diabetic ulcers, cavity wounds, lacerations, skin abrasions, graft and donor sites, post operative surgical wounds, superficial and partial thickness burns. Tegaderm™ Alginate Dressing can be used to control minor bleeding in emergency wounds. Tegaderm™ Alginate Dressing is suitable for use under compression bandaging.			
Contraindications: Do not use to control heavy bleeding. Do not use on patients with a known sensitivity to alginates. Do not for surgical implantation.			
Precautions: Do not use on dry or lightly exuding wounds as the dressing may adhere to the wound. Do not use on full thickness burns. Do not re-use in whole or in part, as it may compromise sterility and/or the performance of the dressing. Tegaderm™ Alginate dressing must not be ingested and must be kept away from children and pets. For use under the avertance of a healthcare professional.			
Directions for Use: Site Preparation and Dressing Application: 1. Debride if necessary and irrigate the wound site in accordance with standard protocols. 2. Remove excess solution from surrounding skin. 3. Select a size of Tegaderm™ Alginate dressing that is slightly larger than the wound. 4. Cut (using sterile scissors) or fold dressing to fit the wound. Apply to wound freely to the dress. 5. Dispose of unused dressing material, according to local regulations and guidelines, due to risk of contamination. For cavity wounds, loosely fill ensuring the dressing does not overlap the wound margins. Cover and secure with a non-occlusive or semi-occlusive secondary dressing, for example Tegaderm™ Hydrocolloid or Tegaderm™ Foam dressing.			
Changing the Dressing and Removal: Tegaderm™ Alginate can remain in-situ for up to 7 days, dependent on patient condition and the level of exudate. 2. Reapply Tegaderm™ Alginate when the secondary dressing has reached its absorbent capacity or whenever good wound care practice dictates. 3. To reapply, gently remove the secondary dressing, if the wound appears dry, saturate the dressing with sterile saline solution prior to removal. 4. Gently remove the dressing from the wound bed and dispose of according to local procedures and guidelines. 5. If required, irrigate the wound site with a suitable wound cleanser prior to application of a new dressing.			
Storage and Shelf Life: This product should be stored at room temperature. Avoid excessive heat and humidity. For shelf life, refer to expiration date which appears on each package.			
How Supplied: Supplied in boxes of individually packaged sterile dressings. If you have any questions or comments, contact your local SA subsidiary. Date of issue: 2020-01			
Symbol	The Title of the Symbol and its Reference Number	The Meaning or Explanatory Text for the Symbol	
ISO Disinfection Numbers (5-17) ISO 15823 – Medical Devices – Symbols Part 1 GB Registered Numbers, batch numbers, lot numbers and lot numbers – Supplied Part			
	5.1.1 MANUFACTURER	Indicates the medical device manufacturer as identified in EU Directive 90/269/EEC, 93/42/EEC and 96/76/EEC.	
	5.1.4 USE BY DATE	Indicates the date after which the medical device is not to be used.	
	5.1.5 BATCH CODE	Indicates the manufacturer's batch code so that the batch and lot can be identified.	
	5.1.6 CATALOGUE NUMBER	Indicates the manufacturer's catalogue code so that the medical device can be identified.	
	5.2.4 STERILIZED USING IRRADIATION	Indicates a medical device that has been sterilized using irradiation (see ISO 13225:5.2.4).	
	5.2.6 DO NOT RESTERILIZE	Indicates a medical device that is not to be re-sterilized.	
	5.2.8 DO NOT USE IF PACKAGE IS DAMAGED	Indicates that if the packaging has been damaged or opened.	
	5.3.4 KEEP DRY	Indicates a medical device that needs to be protected from moisture.	
	5.3.6 PERMIT TAIL OF TEMPERATURE (>25/°C) (<25/°C)	Indicates the temperature limits in which the medical device can be safely applied during a single procedure.	
	5.4.2 DO NOT RE-USE	Indicates a medical device that is intended for one use, or that the device should be discarded after a single procedure.	
	5.4.3 CONSULT INSTRUCTIONS FOR USE	Indicates the need for the user to consult the instructions for use.	
	5.4.4 CAUTION	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for the purposes of this symbol, be presented on the medical device itself.	
	5.1.2 Authorized Representative in European Community	Indicates the authorized representative in the European Community.	
	GREEN DOT TRADEMARK	Indicates a financial packaging recovery company (see European Directive 94/62/EEC, 2002/96/EC, and corresponding national law).	

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