
Mini-Spike®

B | BRAUN



en	Removal and injection spikes for single- and multi-dose vials.	4	no	Aspirasjons- og injeksjonspigger for enkelt- og flerdosebeholdere.	33
de	Entnahme- und Zuspitzspike für Einzel- und Mehrdosenfläschchen.	5	pl	Szpiculce do pobierania i wstrzykiwania do fiolek jedno- i wielodawkowych.	34
bg	Шипове за изтегляне и инжектиране за едно- и многодозови флакони.	7	pt	Remoção e spikes de injeção para frascos de dose única ou multi-dose.	36
cs	Odběrové nebo injekční hroty pro jednodávkové a vícedávkové injekční lahvičky.	8	pt BR	Spikes de remoção e de injeção para frascos de dose única e multidose.	38
da	Fjernelses- og injektionsspyd til hætteglas med en eller flere doser.	10	ro	Vârfuri pentru îndepărtare și injectare destinate flacoanelor unidoză și multidoză.	39
el	Ακίδες αφαίρεσης και έγχυσης για φιαλίδια μίας ή πολλών δόσεων.	12	ru	Наконечники для аспирации и инъекций, предназначенные для однодозных и мультидозных флаконов.	41
es	Punzones de extracción y de inyección para viales de dosis única o múltiple.	13	sk	Aspiračné a injekčné hroty pre jednodávkové a viacdávkové liekovky.	43
et	Ühe- ja mitmeannuseliste vialide eemaldus- ja süsteteravikud.	15	sl	Odstranjevanje in konice za injiciranje za enoodmerne in večodmerne viale.	44
fi	Poisto- ja injektointipiikit yksi- ja moniannosampulleille.	17	sr	Šiljci za ubrizgavanje ili izvlačenje tečnosti iz jednodoznih i višedoznih bočica.	46
fr	Perforateurs de retrait et d'injection pour des flacons mono- et multi-doses.	18	sv	Uttags- och injiceringspetsar för en- och flerdosflaskor.	48
hr	Šiljci za izvlačenje i ubrizgavanje za jednodozne i višedozne bočice.	20	tr	Tek ve çok dozlu flakonlar için çekme ve enjeksiyon spike'ları.	49
hu	Tűskék felszíváshoz/befecskendezéshez egy- és multidózisos ampullákhoz.	21	uk	Аспіраційні канюлі для одноразової та багаторазової аспірації та ін'єкції у флакони.	51
id	Taji atau paku pelepasan dan injeksi untuk vial satu dosis dan multi-dosis.	23	vi	Kim hút và bơm cho các loại ống thuốc một liều hoặc đa liều.	53
it	Dispositivo di prelievo o iniezione per fiale monodose e multidose.	25	zh	用于单剂量和多剂量小药瓶的抽吸和注入穿刺器。	54
kk	Бір және көп дозалы құтыларға арналған алу және енгізу инелері.	26			
lt	Nuėmimo ir injekcijos smaigai vienos dozės ar daugiadoziams flakonams.	28			
lv	Izņemšanas vai injekcijas adatas vienas un vairāku devu pudelītēm.	30			
nl	Verwijderings- en injectiespikes voor single en multidose vials.	31			

en Directions
de Anwendungshinweise
bg Указания за употреба
cs Pokyny pro použití
da Betjeningsvejledning
el Υποδείξεις εφαρμογής
es Aplicaciones
et Kasutusjuhendid

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Device description

Material used:

Housing:	ABS/SAN
Closing Cap:	PP
Protective Cap and Valve Opener:	PE
Filter:	Acrylic copolymer / Nylon support
Valve Disk:	silicone

Mini-Spike® can be used with devices which are in compliance with Luer connectors (e.g., ISO 80369-7).

Mini-Spike® can be used with single- and multi-dose vials and fluid containers of different volumes.

Removal and injection spikes for single- and multi-dose vials.

Sterility

Sterilized by Gamma radiation.

Do not use if package is damaged.

Visually inspect the product packaging to ensure that the sterile barrier system is intact.

Intended purpose

Mini-Spike® is a vented dispensing pin for injecting/withdrawing fluids from single and multi-dose containers (i.e. glass vials, semi rigid and glass IV container).

Do not re-use.

Indication

- Withdrawal and injection from/in multi-dose vials
- Withdrawal and injection from/in multi-dose vials with additional filtration of fluids (Mini-Spike® Filter and some Mini-Spike® Chemo variants only)
- Preparation of ready-for-use injections (according to the SmPC of the drug)
- Preparation of cytostatic injections (according to the SmPC of the drug)

	A	B	C	D	E	*
Mini-Spike® (REF 4550242)	✓					green
Mini-Spike® Micro-Tip (REF 4550510)	✓		✓			green
Mini-Spike® V (REF 4550560)	✓				✓	green
Mini-Spike® Filter (REF 4550234)	✓			✓		blue
Mini-Spike® Filter Micro-Tip (REF 4550528)	✓		✓	✓		blue
Mini-Spike® Filter V (REF 4550579)	✓			✓	✓	blue
Mini-Spike® Chemo (REF 4550340)		✓		✓		red
Mini-Spike® Chemo Micro-Tip (REF 4550536)		✓	✓			red
Mini-Spike® Chemo V (REF 4550587)		✓		✓	✓	red

* Closing cap colour

A: Venting filter 0.45 µm

Protects the medication from contamination by bacteria.

B: Venting filter 0.2 µm

Protects the medication against contamination by bacteria and the user from contamination by aerosols. Specially suited for use with toxic substances.

C: Micro-Tip

Reduces the residual volume.

Specially suited to small volume, expensive medications.

D: Fluid filter 5 µm

Protects against contamination by particles.

Specially suited for use with medications in powder form.

E: Valve

The valve also closes upside down containers and is drip proof.

Patient population

Mini-Spike® is not directly used on patients as it is only used for drug admixture. No gender or age related limitations.

Intended User

Mini-Spike® should only be used by healthcare professionals (e.g., professional nurses, licensed and/or certified physicians, healthcare workers, midwives, paramedics, pharmacists) who have been adequately trained within the educational background in this technique.

Upon medical assessment and adequate instruction, patients and/or caregivers may also be allowed to take over defined handling steps according to national guidelines. Healthcare workers are responsible for ensuring that patients and/or caregivers are instructed in the correct use of Mini-Spike® and that the instruction of caregiver and/or patient is documented by the healthcare worker according to local documentation guidelines.

Contraindications

Withdrawal of cellular containing blood samples.

Reconstitution of solid pharmaceuticals with additional solvent. There are special transfer devices available for this indication.

Residual Risks/Side effects

No side effects known.

The following complications may occur during the procedure of liquid transfer, reconstitution or dilution:

- Microbiological or particulate contamination as a result of poor aseptic technique up to infection
- Medication errors

Precautions

Do not re-sterilize.

Always refer to the SmPCs of the drugs / solutions for admixture procedure.

Warning

Do not re-use. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness, or death of the patient.

Every dispensing pin is intended to be used with a single container and is not made for re-use.

Do not use if packaging or product is damaged or contaminated. Use of damaged or contaminated device might lead to injury, illness, or death of the patient.

The duration of the application is applied in accordance with the used pharmaceutical and its SmPCs as well as in accordance with the national guidelines and/or hospital protocols.

If product is used with cytotoxic or antineoplastic drugs, please mind safe handling instructions and guidelines / institutional protocols.

Do not remove the dispensing pin/spike from the stopper / vial. The injection pin may cause injuries with careless handling.

Solution from a filled syringe into the container must only be injected upright. Vials with inserted Mini-Spike® are to be stored in a vertical position (upright). In case of repeated withdrawal of fluids from a multi-dose vial: disinfect the female luer connection / swabable valve as per institutional protocol.

For variants with Liquid filter / Particle filter: When using the following pharmaceuticals, the absorption is larger than 15% due to the liquid filter. The resulting loss has to be taken into account:

- Corticotropin 30 %
- Cosyntropin 35 %
- Cyanocobalamin 70 %
- Ergonovine Maleate 45 %
- Glucagon 60 %
- Insulin 90 %
- Vasopressin 95 %

Every withdrawal device is intended to be used with a single container and is not made for reuse, because even by applying the device according to its intended use there might occur fissuring or stress cracking due to mechanical and/or chemical influences.

Substances causing stress cracks:

Various inorganic and organic surface-active substances, solvents (e. g. macrogol/polyethylene glycol and/or benzyl alcohol) in medication or active substances (e.g. carboplatin, daunorubicin, fluorouracil, cisplatin, paclitaxel) can lead to stress and tension cracks in plastics. With these aggressive substances we recommend use of the spike for a single aspiration only.

If hygiene rules are disregarded, patient-cross-contamination is potentially possible, in spite of the use of the product.

Operating Instructions

- 1) Remove protective cap from the medication container and disinfect the stopper.
- 2) Remove protective cap from the spike.
- 3) Insert spike.
- 4) Open snap cap.
- 5) Connect syringe and aspirate / inject required volume.
- 6) **Important:** Solution from a filled syringe into the container (Fig. 6) must only be injected right-side up. Risk of excess pressure and impaired deairing.

Duration of use

Please refer to the SmPC of the medication as well as valid hygiene and safety regulations.

Disposal

Disposal acc. to local guidelines and/or clinical protocols.

Storage and handling conditions

Keep away from sunlight.

Keep dry.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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	A	B	C	D	E	*
Mini-Spike® (REF 4550242)	✓					grün
Mini-Spike® Micro-Tip (REF 4550510)	✓		✓			grün
Mini-Spike® V (REF 4550560)	✓				✓	grün
Mini-Spike® Filter (REF 4550234)	✓			✓		blau
Mini-Spike® Filter Micro-Tip (REF 4550528)	✓		✓	✓		blau
Mini-Spike® Filter V (REF 4550579)	✓			✓	✓	blau
Mini-Spike® Chemo (REF 4550340)		✓		✓		rot
Mini-Spike® Chemo Micro-Tip (REF 4550536)		✓	✓			rot
Mini-Spike® Chemo V (REF 4550587)		✓		✓	✓	rot

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Distributed by: **B. Braun Medical Inc.**
B. Braun Melsungen AG