ETHILON™

Polyamide 6 or Polyamide 6,6

Sterile Synthetic Non-Absorbable Surgical Suture

ETHILON™ Polyamide 6 or Polyamide 6,6 Sterile Synthetic Non-Absorbable Surgical Suture	en	ETHILON™ Poliammide 6 o Poliammide 6,6 Sutura chirurgica non assorbibile sintetica sterile	it
ETHILON™ Полиамид 6 или полиамид 6,6 Стерилен синтетичен нерезорбируем хирургичен конец	bg	ETHILON™ Poliamidas 6 arba poliamidas 6,6 Sterilus, sintetinis, nerezorbuojamasis chirurginis siūlas	lt
ETHILON™ Polyamid 6 nebo Polyamid 6,6 Sterilní syntetický nevstřebatelný šicí materiál	cs	ETHILON™ Poliamids 6 vai poliamids 6,6 Sterils, sintētisks, neabsorbējams ķirurģiskais diegs	lv
ETHILON™ Polyamid 6 eller polyamid 6,6 Steril, syntetisk, ikke-resorberbar kirurgisk sutur	da	ETHILON™ Polyamide 6 of polyamide 6,6 Steriel synthetisch niet-resorbeerbaar chirurgisch hechtmateriaal	nl
ETHILON™ Polyamid 6 oder Polyamid 6,6 Steriles, synthetisches, nicht resorbierbares chirurgisches Nahtmaterial	de	ETHILON™ Polyamid 6 eller Polyamid 6,6 Steril syntetisk ikke-absorberbar kirurgisk sutur	no
ETHILON™ Πολυαμίδιο 6 ή Πολυαμίδιο 6,6 Αποστειρωμένο συνθετικό μη απορροφήσιμο χειρουργικό ράμμα	el	ETHILON™ Poliamid 6 lub poliamid 6,6 Jalowe, syntetyczne, niewchłaniałne nici chirurgiczne	pl
ETHILON™ Poliamida 6 o Poliamida 6,6 Sutura quirúrgica estéril sintética no absorbible	es	ETHILON™ Poliamida 6 ou Poliamida 6,6 Fio de sutura cirúrgico esterilizado, sintético e não absorvível	pt
ETHILON™ Polüamiid 6 või polüamiid 6,6 Steriilne sünteetiline mitteresorbeeruv kirurgiline õmblusmaterjal	et	ETHILON™ Poliamidă 6 sau poliamidă 6,6 Fir de sutură chirurgicală steril, sintetic, neresorbabil	ro
ETHILON™ Polyamidi 6 tai polyamidi 6,6 Steriili, synteettinen, resorboitumaton kirurginen ommelaine	fi	ETHILON™ Polyamid 6 alebo Polyamid 6,6 Sterilná syntetická nevstrebateľná chirurgická niť	sk
ETHILON™ Polyamide 6 ou Polyamide 6,6 Fil de suture chirurgicale synthétique stérile non résorbable	fr	ETHILON™ Poliamid 6 ali poliamid 6,6 Sterilna sintetična neresorbilna kirurška nit	sl
ETHILON™ Poliamid 6 ili poliamid 6,6 Sterilni, sintetički, neresorptivni kirurški konac	hr	ETHILON™ Polyamid 6 eller Polyamid 6,6 Steril, syntetisk, icke-resorberbar kirurgisk sutur	sv
ETHILON™ Poliamid 6 vagy poliamid 6,6 Steril, szintetikus, nem felszívódó sebészi varróanyag	hu	ETHILON™ Poliamid 6 veya Poliamid 6,6 Steril Sentetik Emilmeyen Cerrahi Sütür	tr

Instructions for Use

ETHILON™

Polyamide 6 or Polyamide 6,6

Sterile Synthetic Non-Absorbable Surgical Suture

Description

ETHILON™ Suture is a sterile, monofilament, synthetic, non-absorbable surgical suture composed of ≥ 99.0 weight percent polyamide 6 or polyamide 6,6. The empirical molecular formula of the polyamide 6 polymer is $[\text{NH-CO}-(\text{CH}_2)_5]_n$ and the empirical molecular formula of the polyamide 6,6 polymer is $[\text{NH}-(\text{CH}_2)_6-\text{NH-CO}-(\text{CH}_2)_4-\text{CO}]_n$.

ETHILON™ Suture is available undyed and dyed black with up to 1.0 weight percent Hematine HCK (Color Index 75290), and dyed green with up to 0.6 weight percent D&C Green No. 5 (Color Index 61570) to enhance visibility in the surgical field.

ETHILON™ Suture is available in a range of gauge sizes and lengths, non-needled or attached to needles of various types and sizes, and in presentations as described in the How Supplied section.

ETHILON™ Suture complies with the requirements of the European Pharmacopoeia (Ph. Eur.) for Sterile Polyamide 6 or Polyamide 6,6 Suture and the United States Pharmacopeia (USP) for Non-Absorbable Surgical Suture.

The European Pharmacopoeia recognizes units of measure Metric and Ph. Eur. sizes as equivalent which is reflected on the labeling.

ETHILON™ Suture is intended for use only by healthcare professionals who are trained in surgical suturing techniques.

Clinical benefits expected from soft tissue approximation and/or ligation are to promote wound healing by primary intention and to avoid bleeding from vessels during and after surgical procedure.

A summary of safety and clinical performance can be found at the following link (upon activation):

<https://ec.europa.eu/tools/eudamed>.

Indications / Intended Use

ETHILON™ Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurosurgical procedures.

Contraindications

Due to the gradual loss of tensile strength that may occur over prolonged periods *in vivo*, ETHILON™ Suture should not be used where permanent retention of tensile strength is required.

Warnings

Healthcare professionals should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing ETHILON™ Suture for wound closure, as risk of tissue separation / wound dehiscence leading to impaired healing may vary with the site of application and the suture material used.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Like all foreign bodies, ETHILON™ Suture may potentiate infection.

Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation that may result in device failure and/or cross-contamination that may lead to

infection or transmission of blood-borne pathogens to patients and healthcare professionals.

Discard unintentionally opened / partially used / used devices and packages.

Precautions

As with any suture material, adequate knot security requires the standard surgical technique of flat, square ties with additional throws as warranted by surgical circumstance and the experience of the healthcare professional. The use of additional throws may be particularly appropriate when knotting monofilament sutures. Improper knotting or damage to suture during use can cause extended surgery time or treatment failure and additional surgery.

In handling this or any other suture material, care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Healthcare professionals should exercise caution when handling surgical needles to avoid inadvertent needle stick injury that may result in transmission of blood-borne pathogens from contaminated needles. Broken needles may result in extended or additional surgeries or residual foreign bodies. In the event of a product malfunction before use, such as a bent, broken or detached needle, or suture damage, the product should be discarded and a new one obtained to begin the procedure. In the event of product malfunction during use, it is up to the discretion of the healthcare professional whether to continue or discontinue usage of the product and how to complete the procedure.

Adverse Reactions / Undesirable Side Effects

Adverse reactions associated with the use of this device include tissue separation / wound dehiscence leading to impaired healing. Other associated adverse events include gradual loss of tensile strength over time, calculus formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine or bile occurs, minimal inflammatory tissue / foreign body reaction, and transient local irritation at the wound site. Like all foreign bodies, ETHILON™ Suture may potentiate infection.

Suture breakage may result in hemorrhage. Suture material cutting or tearing through the tissue at the time of suture

tensioning may cause soft tissue injury. Improper knotting or damage to suture during use can cause extended surgery time or treatment failure and additional surgery.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

Healthcare professionals should convey adverse reactions, undesirable side effects, and risks associated with the product and the procedure to the patient and advise the patient to contact a healthcare professional in case of any deviation from the normal postoperative course.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

Magnetic Resonance Imaging (MRI) / Carcinogenic, Mutagenic, Toxic to Reproduction (CMR) / Endocrine Disrupting (ED) Safety Information

ETHILON™ Suture is MR safe. No known CMR Category 1a/1b and ED substances are present at >0.1%. Category 1a/1b are defined as known or presumed human carcinogen (H340), mutagen (H350) or reproductive toxicant (H360), based on human evidence and animal studies.

Application / Instructions for Use

Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique, and wound characteristics.

Discard needles in “sharps” containers. Discard unintentionally opened / partially used / used devices and packages according to your facility’s policies and procedures concerning biohazardous materials and waste.

Performance / Actions

ETHILON™ Suture elicits an initial, minimal inflammatory reaction in tissues, that is followed by a gradual encapsulation of the suture by fibrous connective tissue. While polyamide is not absorbed, progressive hydrolysis of the polyamide *in vivo* may result in gradual loss of tensile strength over time.

Sterility

ETHILON™ Suture is sterilized by irradiation. Do not resterilize. Do not use if package is opened or damaged.

Storage

No special storage conditions required. Do not use after expiry date.

How Supplied

Please note that not all sizes are available in all markets. Please contact your local sales representative for size availability.

ETHILON™ Suture is available sterile, monofilament strands in USP sizes 11-0 through 2 (metric sizes 0.1 to 5.0) in a variety of lengths, with and without permanently attached needles.

ETHILON™ Suture is also available in configurations containing the following:

1. Retention tubing (elastomer tubing), that is used to spread the load of the suture at the surface of the skin.

ETHILON™ Suture is available in one, two, or three dozen units per box.

Traceability

The following specific information can be found on the device packaging label: Catalogue number, Batch Code, expiry and manufacturing date, manufacturer name, address and the website and a Unique Device Identification bar code with the Global Trade Item Number information.

Symbols Used on Labeling

Catalogue number

Medical Device

Do not use if package is damaged and consult instructions for use

Caution

Do not re-use

Do not re-sterilize

Single sterile barrier system with protective packaging inside

Sterilized using irradiation

Batch code

Date of manufacture

Use-by date

Unique Device Identifier

Manufacturer

Authorized representative in the European Community/European Union

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Packaging unit

PEEL/LIFT/OPEN HERE

TEAR HERE

www.e-ifu.com
EU: +800 8888 2020
EU: +32 2 4037222

Consult instructions for use or consult electronic instructions for use.

EU: Call paper on demand help desk to get paper copies free of charge within 7 days.

MR Safe

Dyed - Non-absorbable - Monofilament - Uncoated - Suture

Undyed - Non-absorbable - Monofilament - Uncoated - Suture

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AB: Basılı kopyaları 7 gün içinde ücretsiz olarak talep üzerine almak için danışma hattını arayın.
Ethicon, LLC
Johnson & Johnson Medical GmbH