# Gazin® RK

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## en Instructions for Use

**Product description**

Gazin® gauze balls RK, non-sterile, consist of 20-thread count cotton gauze compliant with EN 14079. The X-ray detectable thread contained in the gauze balls consists of barium sulphate filaments and is wound with polypropylene.

**Product composition**

100% cotton, X-ray detectable thread (blue, polypropylene, polyester thread, barium sulphate).

**Intended purpose**

Gazin® gauze balls are used for absorbing body fluids during surgical procedures.

**Usage information/Precautions**

This product is non-sterile. Gazin® gauze balls must be sterilised prior to surgically invasive use; resterilisation is not permitted. The product is intended for single use only.

**Sterilisation instructions**

**1. Device requirements**

Only sterilisers, which comply to the requirements of the norms EN 285 or EN 13060 and which use a pre-vacuum system for the sterilisation are allowed to be used.

The use of the gravitational method is not allowed in the EU.

**2. Sterilisation preparation**

Before the product is sterilised, a suitable sterile barrier system must be selected, validated and used by the person responsible for the sterilisation (sterilising technician) in accordance with ISO 11607. The sterile barrier must be big enough for the product so that the seal is not under tension. The retail packaging is unsuitable for sterilisation.

**3. Sterilisation**

The sterilisation procedure for the product must be validated by means of moist heat at 134°C for 7 minutes in accordance with ISO 17665. Do not use after sterilisation if package is damaged or accidentally opened.

**4. Storage**

After removing the required quantity of Gazin® gauze balls, non-sterile, seal the packaging and store dry and protected from dust. Sterilised products must be stored dry and protected from dust in microbe-resistant sterile packaging.

**Additional information**

The above instructions were assessed by Lohmann & Rauscher as suitable for preparation for sterilisation. The sterilising technician is responsible for ensuring that the sterilisation carried out with the employed equipment, materials and staff achieves the required results. To this end, validation of the sterilisation process and routine monitoring of the procedure is normally required. Deviations from the procedure are the responsibility of the sterilising technician. National requirements for steam sterilisation must be considered by the user.

**Disposal**

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

**General instructions**

In case of a serious incident, contact the manufacturer and the responsible health authority.

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