Instructions for Use

TWINFIX[◇] Ti Preloaded Suture Anchors

EN English

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°Trademark of Smith & Nephew, registered U.S. Patent & Trademark Office. This product may be covered by one or more U.S. patents. See smith-nephew.com/patents for details.
©2006, 2024 Smith & Nephew, Inc. All rights reserved.

06/2024 10600147 Rev. J

TWINFIX[◊] Ti Preloaded Suture Anchors

Instructions for Use

Device Description

Smith & Nephew TWINFIX[◊] Ti Preloaded Suture Anchors consist of a suture anchor with attached non-absorbable suture(s) assembled to an insertion device. The suture anchor is intended to provide secure reattachment of soft tissue to bone. Attachment of soft tissue is performed by the surgeon's preferred technique.

Performance Characteristics

The TWINFIX Ti Preloaded Suture Anchors are 2.8 mm, 3.5 mm, 5.0 mm, and 6.5 mm diameter anchors with a helical screw design for providing proper suture fixation. The anchors are supported by an internal driver for delivery and insertion.

How Supplied

The TWINFIX Ti Preloaded Suture Anchors are supplied sterile, for single use only.

Contents

- 1 ea Suture anchor
- 1 ea Insertion device
- 1 or 2 ea Non-absorbable suture(s)

Refer to individual suture anchor product label(s) for suture size and quantities.

Materials

- **Suture anchor** – 100% titanium alloy
- **Insertion device** – Stainless steel shaft with polycarbonate and ABS handle
- **Suture** – Braided, uncoated, UHMW polyethylene and UHMW polyethylene with monofilament polypropylene cobraid, non-absorbable

Intended Users

TWINFIX Ti Preloaded Suture Anchors are intended for use by healthcare professionals in accordance with these instructions for use. The use environment is a professional healthcare facility.

Intended Use

TWINFIX Ti Preloaded Suture Anchors are intended for use only for the reattachment of soft tissue to bone.

Indications for Use

TWINFIX Ti Preloaded Suture Anchors are indicated for use only for the reattachment of soft tissue to bone for the following indications:

Shoulder

- Bankart lesion repairs
- SLAP lesion repairs
- Rotator cuff tear repairs
- Capsular shift or capsulolabral reconstructions
- Biceps tenodesis

Knee

- Extra-capsular repairs:
 - medial collateral ligament
 - lateral collateral ligament
 - posterior oblique ligament

Contraindications

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Surgical procedures other than those listed in the Indications for Use section.
- Pathological conditions of bone, such as cystic changes or severe osteopenia, which would compromise secure anchor fixation.
- Pathological conditions in the soft tissues to be attached that would impair secure fixation by suture.
- Comminuted bone surface, which would compromise secure anchor fixation.
- Physical conditions, which would eliminate, or tend to eliminate, adequate anchor support or retard healing.

Target Patient Group

The patient population is comprised of adult patients for whom the reattachment of soft tissue to bone during orthopedic procedures is determined to be necessary by the prescribing healthcare professional.

Expected Clinical Benefit

TWINFIX Ti Preloaded Suture Anchors facilitate successful reattachment of soft tissue to bone as demonstrated through improvement in postoperative scores.

Links to the Summary of Safety and Clinical Performance (SSCP)

For use in the European Union (EU), a summary of safety and clinical performance (SSCP) has been developed in conformance with the EU regulations. Specific information about the subject devices, including references/catalog numbers, manufacture name, and trade name, may be used to search for the appropriate SSCP on the European database on medical devices (EUDAMED). Multiple elements of information may be needed to adequately identify the device in EUDAMED. This is the SSCP location after the launch of EUDAMED: <https://ec.europa.eu/tools/eudamed>

Warnings

- **Do not use if the product sterilization barrier or its packaging is compromised.**
- **TWINFIX[®] Ti Preloaded Suture Anchors are supplied sterile and are for SINGLE USE ONLY. Components are provided in sterile packaging and contents will remain sterile unless package is opened or damaged. Discard any open, unused products and do not use after expiration date. Do not clean, resterilize, or reuse the device, as this may damage or compromise the performance resulting in product malfunction, failure, or patient injury. Cleaning, resterilization, or reuse of the implant may also expose the patient to the risk of transmitting infectious diseases.**
- **It is the healthcare professional's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.**
- **Read these instructions completely prior to use.**
- **Product must be stored in the original sealed pouch.**
- **Titanium alloys contain elements that may stimulate allergic hypersensitive responses by the immune system. These elements are titanium, aluminum, and vanadium (Ti, Al, and V). When sensitivity is anticipated, appropriate preoperative testing should be conducted.**
- **Incomplete anchor insertion may result in poor anchor performance.**
- **Breakage of suture anchor can occur if predrilling is not performed prior to implantation.**
- **Do not attempt to implant this device within cartilage epiphyseal growth plates or non-osseous tissue.**
- **Associated instruments for TWINFIX Ti Preloaded Suture Anchors are sold separately and are provided NON-STERILE. These instruments must be properly cleaned and sterilized prior to use.**
- **Inspect the drill bit for damage prior to use. Replace a damaged or worn drill bit. Do not attempt to straighten or sharpen. Sharpening the drill bit will alter the implantation site and could affect anchor stability.**
- **Do not resterilize or reuse anchors, sutures, insertion devices.**

- **As with any surgical procedure, there is risk involved. Potential complications accompanying such implant surgery may include, but are not limited to, the following:**
 - **Breakage of the suture**
 - **Loss of fixation or pullout of suture anchors**
 - **Mild inflammatory reaction**
 - **Foreign body reaction**
 - **Infection, both deep and superficial**
 - **Allergic reaction**

Precautions

R_x U.S. Federal law restricts this device to sale by or on the order of a physician.

- Hazards associated with reuse of this device include, but are not limited to, patient infection and/or device malfunction.
- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- Do not use sharp instruments to manage or control the suture.
- Hard bone conditions require preparation by predrilling the insertion site to reduce the potential of torsional overload. Predrilling extracts the core diameter of the suture anchor and creates a countersink broach for insertion of the device tip. Predrilling with the appropriate drill bit is the preferred method of site preparation.
- Ensure the anchor placement is aligned with the drilled hole. Proper alignment is essential for successful repair.
- Excessive force during insertion can cause failure of the suture anchor or insertion device. A two-finger AO technique should be used to insert the anchor.
- As in all suture techniques, the suture is intended to approximate soft tissue to bone for a length of time appropriate for biological attachment of soft tissue to bone. The suture is not intended to provide indefinite biomechanical integrity.
- Bone must be adequate to allow proper placement of suture anchor.
- Users shall not alter implant or instrumentation, otherwise performance may be compromised.
- Postoperative range of motion is to be determined by the physician.
- After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.

Magnetic Resonance Imaging (MRI) Safety

The TWINFIX[°] Ti Preloaded Suture Anchors are manufactured from titanium alloy. Nonclinical testing has demonstrated that the suture anchors are MR Conditional. MR Conditional refers to an item that has been demonstrated to pose no known hazards in a specified MR environment under specified conditions of use. This implant can be scanned safely under the following conditions:

- Static magnetic field of 1.5 Tesla (1.5 T) or 3 Tesla (3.0 T)
- Maximum spatial field gradient of 2,000 G/cm (20.0 T/m)
- When scanning in a 1.5 T or 3 T scanner, patients can be scanned with a whole body averaged SAR of 2.0 W/kg (Normal Operating Mode) without landmark restrictions.
- The scanner SAR restrictions above apply to whole-body RF coil using a circularly polarized transmit mode. No other RF coil types (e.g. extremity, head, neck) were tested.

1.5 T RF Heating

In non-clinical testing under scan conditions defined above, the suture anchors are expected to produce a temperature rise less than or equal to 3.1° C at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg, as assessed by calorimetry for 15 minutes of continuous scanning in a 1.5 T Siemens Espree (MRC30732) MR scanner with SYNGO MR B19 software.

3 T RF Heating

In non-clinical testing under scan conditions defined above, the suture anchors are expected to produce a temperature rise less than or equal to 4.0° C at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg, as assessed by calorimetry for 15 minutes of continuous scanning in a 3 T Siemens Trio (MRC20587) MR scanner with SYNGO MR A35 4VA35A software.

MR Artifact

In non-clinical testing, the image artifact caused by the suture anchor extends radially approximately 1.8 cm from the device when scanned in nonclinical testing using a gradient echo sequence in a 3 T Siemens Trio (MRC20587) MR scanner with SYNGO MR A35 4VA35A software and body coil.

Adverse Events Reporting

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

complaints@smith-nephew.com

Instructions for Use

1. Using the appropriate Smith & Nephew drill bit (sold separately), prepare the anchor insertion site (refer to the Use of Smith & Nephew Drill Bits section of this document). Either predrilling or use of a punch-type device are recommended for site preparation, depending on the bone quality. Predrilling is recommended if the quality of the bone is unknown.

CAUTION: Hard bone conditions require preparation by predrilling the insertion site to reduce the potential of torsional overload. Predrilling extracts the core diameter of the suture anchor and creates a countersink broach for insertion of the device tip. Predrilling with the appropriate drill bit is the preferred method of site preparation.

2. Using routine sterile technique, remove the TWINFIX[°] Ti Preloaded Suture Anchor with attached insertion device from the packaging.
3. Establish axial alignment of the suture anchor to the drill hole. While stabilizing the insertion site, rotate the anchor clockwise into the hole with the insertion device using a two-finger AO technique. Continue rotating the anchor until desired placement is achieved by visual inspection.

Notes:

- The first horizontal indicator line on the tip of the insertion device identifies the top of the anchor.
- The second horizontal indicator line is a 1 cm depth reference mark.
- Axial lines on the insertion device indicate suture orientation within the anchor.

CAUTION: Excessive force during insertion can cause failure of the suture anchor or insertion device. A two-finger AO technique should be used to insert the anchor.

If more torque is required to insert the anchor, **STOP** and ensure that the anchor size, drill hole size, and depth are correct for the bone conditions encountered. It may be necessary to reduce the anchor size or increase the drill hole size to achieve optimal insertion force. It is the responsibility of the surgeon to determine the patient's bone condition, appropriately prepare the insertion site, and determine the suitability of the implant for the procedure.

4. Disengage the suture anchor from the inserter:
 - a. Hold the handle securely in the palm and pull back on the finger grips using two fingers to retract the suture tensioning mechanism.
 - b. While holding the tensioning mechanism in the back position with the finger grips, slowly remove the inserter from the insertion site until the suture tips are visible. The suture will release and feed through the inserter as it is removed.

Alternatively:

- a. Pull back on the finger grips and unwind the suture with the free hand.
- b. Slowly remove the inserter.
- c. As the inserter is removed, pinch the suture below the inserter using the free hand to verify suture release from the inserter.

CAUTION: Do not use sharp instruments to manage or control the suture.

5. Discard the insertion device.

Anchor Removal

If it becomes necessary to remove the anchor from the implant site, reengage the anchor head completely with the insertion device and rotate the anchor in a counterclockwise direction. If the suture is still attached, it can be pulled through the insertion device shaft using a suture passer. This will provide counter-force to the anchor.

Instructions for Reassembly

If circumstances necessitate reassembling the suture and anchor to the insertion device, such as surgeon preference for alternate suture with the anchor, or partial deployment of the anchor prior to insertion, reassemble as follows:

1. Disassemble the suture anchor from the insertion device:
 - a. While retracting the suture tensioning mechanism using the finger grips, gently pull the anchor and suture from the inserter shaft.
2. Remove the suture from the anchor:
 - a. Gently pull the suture through the anchor eyelets.

Note: Suture removal is only necessary if an alternate suture is required.
3. Affix the desired alternate suture to the anchor:
 - a. Feed 36 in. (91 cm) minimum length of suture(s) through the anchor eyelet(s).
 - b. Equalize the length of each suture leg.
4. Assemble suture anchor into the inserter:
 - a. Feed the suture(s) into the tip of the inserter shaft using a 12 in. (30 cm) suture passer until the suture tips extend through the handle.
 - b. Pull the suture through the handle until the anchor is near the tip of the inserter shaft.
 - c. Align the anchor eyelets to the axial lines on the tip of the inserter shaft and seat the anchor into the tip.
5. Set the suture tension:
 - a. Wrap the suture around the suture posts on the handle until the length of the suture tail is equal to the handle length plus 1/2 in. (1.3 cm).
 - b. Retract the suture tensioning mechanism using the finger grips.
 - c. Slide the suture tails into the suture groove on the handle.
 - d. Release the finger grips to set the suture tension.

The device is now ready to use.

Use of Smith & Nephew Drill Bits (sold separately)

The following drill bit size is recommended for the following anchors:

| Drill Bit | Anchor |
|-----------|-------------------|
| 1.8 mm | 2.8 mm |
| 2.5 mm | 3.5 mm and 5.0 mm |
| 3.5 mm | 6.5 mm |

- A preferred means of drilling should be selected by the surgeon, i.e., power drill, manual drill, etc.
- Proper hole depth is achieved when the shoulder on the proximal end of the drill bit contacts the bone surface and drilling cannot continue.

WARNING: Inspect the drill bit for damage prior to use. Replace a damaged or worn drill bit. Do not attempt to straighten or sharpen. Sharpening the drill bit will alter the implantation site and could affect anchor stability.

Lifetime of Device

The functional lifetime of the TWINFIX[®] Ti Preloaded Suture Anchors is 12 weeks, which takes into consideration the primary healing period for the reattachment of soft tissue to bone.

Safe Disposal

After use, this device may be a potential biohazard/sharps hazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.

Compatible Devices

| REF no. | Drill Bit | Description |
|---------|-----------|---|
| 7209506 | 1.8 mm | Drill for TWINFIX Ti 2.8 mm anchor |
| 7209502 | 2.5 mm | Drill for TWINFIX Ti 3.5 mm and 5.0 mm anchor |
| 7210215 | 3.5 mm | Drill for TWINFIX Ti 6.5 mm anchor |

For Further Information

If further information, including warranty, on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative.

Symbols Glossary

| | | | | | |
|---|--------------------------|---|---|---|---|
|  | Manufacturer |  | CE Mark |  | MR conditional |
|  | Date of Manufacture |  | Do not use if package is damaged. Consult Instructions for Use. |  | Consult electronic instructions for use |
|  | Use-by date |  | Sterile, sterilized using ethylene oxide |  | Prescription use only |
|  | Catalog Number |  | Do not re-use |  | Medical Device |
|  | Batch Code |  | Do not re-sterilize |  | Authorized Representative in the European Community |
|  | Unique Device Identifier |  | Single sterile barrier with protective packaging inside |  | European Importer |
|  | Caution | | | | |

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Smith & Nephew Operations B.V.
©2006, 2024 Smith & Nephew, Inc. Tüm hakları saklıdır.