CRPslide

A RAPID SLIDE TEST FOR THE QUALITATIVE UND QUANTITATIVE DETERMINATION OF C-REACTIVE PROTEIN IN HUMAN-SERUM

Principle:

The CRP reagent is a suspension of polystyrene latex particles coated with the gamma globulin fraction of antihuman CRP specific serum. When CRP is present in the sample, presence of agglutination indicates a content of CRP equal or greater than 6 mg/l, without previous sample dilution.

Summary:

C-reactive protein is the classic acute phase protein to inflammatory reactions. It is synthesized by the liver and consists of five identical polypeptide chains forming a five-membered ring of molar mass 120.000 daltons. CRP is the most sensitive of the acute phase reactants, and its concentration increases rapidly during inflammatory processes. Complexed CRP activates the complement system beginning with C1q. CRP then initiates opsonization and phagocytosis of invading cells, but its main function is to bind and detoxify endogenous toxic substances produced as a result of tissue damage. CRP assays are used to detect systematic inflammatory processes (apart from certain types of inflammation such as SLE and Colitis ulcerosa); to assess treatment of bacterial infections with antibiotics; to detect intrauterine infections with concomitant premature amniorrhexis; to differentiate between active and inactive forms of disease with concurrent infection, e.g. in patients suffering from SLE or Colitis ulcerosa; to therapeutically monitor rheumatic diseases and assess anti-inflammatory therapy; to determine the presence of postoperative complications at an early stage, such as infected wounds, thrombosis and pneumonia; and to distinguish between infection and bone marrow transplant rejection. A variety of methods, such as nephelometry and turbidimetry are available for the determination of CRP. This CRP assay is based on the principle of immunological agglutination.

Reagent Concentration:

R1:

Latex particles coated with anti-CRP preservative

R2:

Protein containing solution with C-reactive protein > 6 mg/l preservative

R3:

Protein containing solution with C-reactive protein < 6 mg/l preservative

Preparation and stability

The reagents and control sera are stable:

Up to the expiration date at +2°C to +8°C

Do not freeze!

Specimen:

Collect serum using standard sampling tubes

Use fresh clear serum specimens. Plasma, lipemic serum or microbial contamination may cause erroneous results. If the test cannot be performed immediately, store the specimen at +2 to +8°C for up to 24 h. For longer storage, freeze the serum.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents. Contaminated sera and a longer reaction time than 3 min. may cause false positive agglutination.

Precautions:

Reagents containing sodium acid may combine with copper and lead plumbing to form highly explosive metal azides. Dispose of reagent by flushing with large amounts of water to prevent azide build up. The positive and negative controls were prepared from human sera which have been tested using FDA licensed methods and found to be non reactive for HbsAg and HIV antibodies. However, no test method can offer complete assurance that infectious agents are absent. Therefore all human specimens should be considered potentially infectious.

Procedure Qualitative Test:

1. Bring reagents and specimens to room temperature before use.
2. Resuspend the Latex Reagent gently.
3. Place one drop (40 µl) of each test specimen or as reference one drop (40 µl) of the pos. / neg. control on to the reaction slide.
4. Add one drop Latex Reagent (40 µl) to each test field. Use spatula dispenser to spread reaction mixture over entire test field.
5. Rotate the slide for 2 minutes and read immediately under direct light.

Procedure Quantitative Test:

1. Bring reagents and specimens to room temperature before use.
2. Resuspend the Latex Reagent gently.
3. Using saline, dilute the specimens 1:2, 1:4, 1:8, 1:16 or as needed.
4. Place one drop (40 µl) of each test specimen or as reference one drop (40 µl) of the pos. / neg. control on to the reaction slide.
5. Add one drop of the Latex Reagent (40 µl) to each test field. Use stick to spread reaction mixture over entire test field.
6. Rotate the slide for 2 minutes and read immediately under direct light.

Results Qualitative Test:

A negative reaction is indicated by a uniform milky suspension with no agglutination as observed with the negative control (2). A positive reaction is indicated by any observable agglutination in the reaction mixture. The specimen reaction should be compared with the positive control (1).

Results Quantitative Test:

The titer of the serum is the reciprocal of the highest dilution which exhibits a positive reaction.

Interpretation:

positive (1)

Agglutination
within 2 minutes

negative (2)

no Agglutination
within 2 minutes

Reference values:

Adults: < 6 mg/l

CRP is a protein which serum levels increase significantly in an acute phase status as it happens in the most forms of different tissue injuries. It was described the presence of such protein in normal persons at a titer not exceeding 6 mg/l. However the increase occurs in a non specific way in different aggressions like malignant tumor, rheumatic fever, myocardial infarct, inflammations, etc.

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range. For diagnostic purposes the CRP results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

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