# **visocor**® **om**50

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Table of contents

Technical information	22
. Failure and error messages	22
. Customer service	24
. Technical data	25
. Original spare parts and accessories	27
. Technical inspection / calibration check	27
. Explanation of symbols	28
7. <u>Disposal</u>	28
Maintaining the unit	29
Warranty	30

#### Intended use

Thank you for choosing the visocor® OM50 upper arm blood pressure monitor (also referred to in the following as the unit).

It's determined for measuring blood pressure and pulse at the upper arm, those of people 12 years of age and over. This unit is recommended for patients with unstable blood pressure for domestic self administration blood pressure control and treatment support.

The unit has a build in electronic for collection, preparation, analysis and supervision of pulse waves (oscillometric measurement).

Using the pressure variations caused by the pulse at falling cuff pressure, the micro-computer calculates the values systole, diastole, and pulse. Afterwards, the cuff is deflated completely.

A measurement value memory stores the most recent 60 measurement results in each case.

These instructions are intended to help you use the unit reliably and effectively. Keep the instructions together with the unit at all times. The unit must be used in accordance with the methods detailed in these instructions and may not be used for other purposes. Read these instructions carefully before using the unit.

#### 1. Important instructions for patients

Blood pressure measurements on children require specialist knowledge! Please consult your doctor
EN-34

if you wish to measure the blood pressure of a child. Do not under any circumstances use the unit on a baby/infant.

- The cuff has been specially developed for this unit and must not be used with other units.
- The measurement results of automatically measuring blood pressure monitors may be falsified by pregnancy, irregular heartbeat or arteriosclerosis. Measure your own blood pressure in consultation with your doctor.
- Do not under any circumstances place the cuff on or over any critical point, e.g. wound, aneurysm, etc. Risk of injuries! The supply via an intravascular access (infusion) or other medical monitoring devices might be interrupted.
- Do not use the unit without first consulting your doctor if you are undergoing dialysis treatment or taking anticoagulants, platelet aggregation inhibitors or steroids. Internal bleeding may be caused under these conditions.
- In rare cases the cuff can cause bruises on the skin.
- Make sure the air hose is not kinked. A kinked air hose can prevent the cuff from deflating and thereby interrupt the blood flow in the arm for too long.
- The unit contains small parts which could be swallowed by children. It should therefore be kept out of the reach of children at all times. The airhose poses a risk of strangulation. Please follow instructions regarding the handling of cuff to avoid such risk.

## **B** Safety instructions


- If you have had a mastectomy, do not measure blood pressure on the side of the body where the breast has been removed together with the lymph nodes from the armpit.
- Measuring your own blood pressure does not constitute treatment. On no account modify of your own accord the dosage of drugs/medication prescribed by your doctor.
- Please refer to the chapter "Important instructions for use" (Page EN-40) before conducting any measurements yourself.

#### 2. Important technical details

- A consistently good power supply to the unit is essential for trouble-free blood pressure measurements.
- Please use long-life alkaline batteries only (LR6).
- When replacing batteries, always replace all the batteries at the same time.
- You will need 4 x 1.5 V AA/LR6 batteries. Rechargeable batteries have a voltage of only 1.2 V and are thus not suitable
- When using the unit with a mains adapter, please use only the type A1 mains adapter specially tested for medical devices.
- Make sure when operating the unit with the mains adapter that the unit can be disconnected from the mains power supply at all times.
- If the unit is not used for an extended period of time, please remove the batteries. Essentially, all batteries can leak.
- In adverse circumstances, such as in the case of a battery short circuit, the battery compartment cover may heat up.


• If the unit has been stored at temperatures below 0°C, let it stand at room temperature for at least 1

- hour before use.
- The unit must only be operated with original parts. The warranty will be invalidated if the unit is damaged by non-approved accessories!
- Using the unit in the immediate vicinity of mobile phones, microwave appliances or other devices
  with strong electromagnetic fields may result in impaired functioning. Maintain a minimum distance
  of 3 m from such devices when using this unit.
- The display of the pulse frequency is not suitable for checking the frequency of cardiac pacemakers. Cardiac pacemakers and blood pressure monitors do not influence each other in their mode of operation.
- Never open or modify the unit or the cuff this is a medical device and may only be opened by authorised specialists. (Exception: replacing batteries). If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution.
- The cuff may be inflated on the arm only.
- Please comply with the specified operating conditions for measurement. See Technical data, Page EN-54.
- The inflating and measuring procedure can be stopped by pressing the Start/Stop button. In this case, the unit stops the inflation procedure and deflates the cuff.

## Operating the unit

#### 2. Displays

- 1 Memory location or date/time
- 2 SYS = Systole (upper blood pressure value)
- 3 DIA = Diastole (lower blood pressure value)
- 4 PUL 1/min = Pulse calculated pulse frequency per minute
- Pulse signal display or irregular pulse waves (Page EN-47)
- 6 Battery check display
- 7 Deflating cuff
- 8 Inflating cuff
- 9 Memory code

Failure and error messages (Page EN-51)

EN-38

#### 3. Important instructions for use

- Refrain from drinking alcoholic or caffeinated beverages and smoking at least one hour before measuring.
- Please rest at least 5 minutes before measuring. Depending on the previous degree of stress/exertion, this can even require up to one hour.
- Expose the upper arm. Clothing must not impair the blood flow in or from the arm, as this influences the blood pressure at the measurement point.
- Body posture must be relaxed:
  - For this purpose, sit at a table (if possible the height of a dining table, not a coffee table!).
    Rest your back against the backrest of the chair.
  - Rest your forearm completely flat, with the palm of your hand facing upwards.
  - Place your feet on the floor and do not cross your legs.
- It is absolutely essential to sit still and remain quiet. Irregular movements, vibrations (shaking), talking and breathing heavily will affect and may even falsify the measurement. Look out for irregular pulse signals in the display; if necessary, repeat the measurement under more favourable conditions.
- Unusual measurement results are possible with all automatically measuring blood pressure measurement systems from time to time. Check yourself: did you comply with the instructions for use above?
   If necessary, repeat the measurement after allowing the blood circulation in the arm to recuperate briefly; rest for about 3-5 minutes for this purpose.

Operating the unit

Tip: Remain seated, your monitor will switch off automatically about 3 minutes after a measurement.

Afterwards, we recommend that you repeat the measurement.

- In order to rule out differences between left and right and to obtain comparable measurement results, it is important always to take blood pressure measurements on the same arm. Ask your doctor which arm he/she recommends you take your measurements on.
- Blood pressure is not a fixed value. It may fluctuate upwards or downwards by more than 20 mmHg on patients within a few minutes.

#### 4. Initial operation of the unit

Insert the supplied batteries into the unit, if these have not been inserted already.

If the unit is to be operated from the mains power supply, the cable plug of the mains adapter (not included in the scope of delivery) must be inserted into the connection socket of the unit. The batteries are switched off automatically. Please use only the UEBE type A1 mains adapter. See original spare parts and accessories, Page EN-56.

#### 5. Inserting/replacing batteries

Opening the battery compartment:
 Open the cover of the battery compartment on the underside of the unit.

- Inserting batteries:
   Remove the old batteries from the unit and insert the new ones. Please observe proper polarity (marking in battery compartment).
- Closing the battery compartment:
   Close the battery compartment by clipping the battery cover back into the unit.

When the display shows the "Flat battery —" symbol for the first time, the unit will still have enough charge for about 30 measurements. Please change the batteries during this period.

#### 6. Activating and setting the date/time

The unit has a date/time function which you can activate where necessary. This function is deactivated when the unit is delivered.

#### Activating the date/time function

To activate the date/time function, press the Start/Stop button and then immediately, while the full

display is visible, the Memory button. This function can be deactivated at any time in the same way.

#### Setting the date/time


The four-digit number representing the year appears in the display first. Press the Memory button to set the year. Press the Start/Stop button to confirm the selected year and proceed to setting the date.

Set the month with the Memory button and confirm this with the Start/Stop button. Set the day, hours and minutes in the same way.

#### Changing the date/time

To change the date/time function, take out a battery and wait until the display goes blank. Now you can reactivate the date/time function and reset the date and time again.

#### Deactivating the date/time function

To deactivate the date/time function, press the Start/Stop button and then immediately the Memory button.

When the date/time function is activated, the time is also shown in the display when the unit is switched off

#### 7. Attaching the cuff


- Expose the upper arm.
- Push the cuff onto the arm (Fig. 1) until the lower edge of the cuff is 2-3 cm above the crook of the arm (Fig. 2).
- When using the unit on the left arm, the air hose runs to the unit through the middle of the crook of the arm so that the 4 cm long artery marking is automatically located centrally above the pulse position (Fig. 3).


Fia. 1


Fig. 2

Operating the unit

- Bend your arm slightly.
- The upper arm muscle is tensioned slightly.
- This slightly increases the circumference of the arm.
- Now pull the free cuff end tight and close the Velcro fastener (Fig. 4)
- Check if the measurement arrow on the cuff is inside the "marking for arm circumference" at the edge
  of the cuff.
- Lay your extended and relaxed arm with the cuff on the dining table and do not move it during the measurement; do not talk. The palm should be facing upwards.
- Connect the cuff connector to the cuff socket on the left-hand side of the unit. Make sure the connector is fully inserted into the unit.


EN-45

Fig. 4

EN-44

#### 8. Measuring blood pressure

- Switch on the unit by pressing the Start/Stop button. The following display is visible.
- After the unit has completed its calibration against the ambient air pressure, the cuff is automatically inflated to approx. 180 mmHq. If that cuff pressure is not sufficient for measuring the blood pressure level, the device automatically increases the cuff pressure in steps of 30 mmHq until a sufficient cuff pressure has been reached.

- For an inflation pressure exceeding 180 mmHg the inflation process can also be controlled manually: Keep the Start/Stop button pressed immediately after the inflation has started and release it when the cuff pressure has reached 40 mmHg above the expected maximum systolic blood pressure. The inflation stops with releasing the button.
- The actual measurement takes place with the start of deflation. The "♥" symbol in the display flashes until the pulse frequency is displayed.

#### Operating the unit

- The measurement being completed is indicated by a long beep. The determined values for systole, diastole, and pulse are shown on the display.
- The unit switches off automatically after about 3 minutes.

#### 9. Irregular pulse waves and cardiac arrhythmia

If the "w" symbol for irregular pulse waves flashes after a measurement, the unit has recorded irregular pulses during the measurement. This can be set off by an irregular heartbeat (arrhythmia), disturbance caused by movement, talking or even breathing in deeply. The symbol is stored with the respective measurement.

If the symbol appears more frequently, this must be discussed with your doctor! Measurement with irregular pulse waves should therefore be viewed as critical and repeated under more favourable conditions.

EN-46 EN-47

#### 10.Using the memory

Measured results are automatically stored in the memory. The memory can store up to 60 results and the average value. When more than 60 measured values have been stored, the oldest value (No. 60) is deleted to allow the latest value (No. 1) to be recorded.

#### Recalling data


Press the Memory button to call up data. The average value of the stored results is displayed with "A".

When the Memory button is pressed again, the value measured last (No. 1) appears.

Press the Memory button repeatedly to call up further measured values.

When the date/time function is activated, the measured values are displayed alternately with memory number, date and time.

The stored data remain visible in the display for around 30 seconds. The unit then switches itself off.

## Operating the unit

#### Deleting data

To delete individual values, press the Memory button several times to obtain the measured value required.

Then press the Memory button again and hold it down for 8-10 seconds until the value first flashes and then disappears.

To delete the entiry memory, recall the relevant average value display (A) and press the Memory button again, holding it down until the mean value first flashes and then disappears.

EN-48

## What you should know about blood pressure

#### 1. Systolic and diastolic blood pressure values

The cardiovascular system has the important function of supplying all organs and tissues in the body with sufficient amounts of blood and of transporting metabolites. To do so, the heart contracts and expands at a regular rate of about 60 to 80 times per minute. The pressure of the flowing blood on the artery walls caused by the heart contracting is termed systolic. The pressure in the ensuing relaxation phase, when the heart refills with blood, is termed diastolic. During daily measurement you determine both values.

#### 2. Why you measure different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. It can be affected by even slight changes. This explains why values measured by a doctor or pharmacist are higher than those measured at home in the environment you are used to. Changes in the weather, climate changes, or physical or psychological stress can have effects as well.

#### 3. Why you should measure blood pressure regularly

Even the time of day has an influence on your blood pressure. During the day the values are generally higher than during the periods of rest at night. One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is possible only when measurements are taken regularly. Discuss the measurement values with your doctor.

## What you should know about blood pressure

## D

#### 4. What are normal values for blood pressure?

The World Health Organisation (WHO) established the following limits in mmHg (millimetre hydrargy-rum) to classify the blood pressure values.

WHO 2003	Systolic pressure = Upper value mmHg**	Diastolic pressure = Lower value mmHg**
High pressure	from 140*	from 90*
Limit	120 to 139	80 to 89
Normal value	below 120	below 80

<sup>\*</sup> It is sufficient for one of the values to be raised. \*\* millimetres of mercury

EN-50

## Technical information

#### 1. Failure and error messages

Failure encountered	Possible cause	Corrective action
Display Err - 300	Excessive cuff pressure.	- Repeat measurement
	Cuff inflated to maximum due to mo-	- Do not move arm
	vement of the arm or body.	- Do not talk
	The hose may be kinked or the air	Check that hose with air connector is clear
	connector blocked.	(unobstructed)
Display Err - 1	Measurement error. Measurement	- Repeat measurement
	could not be carried out due to	- Do not move arm
	movement of the arm or body.	- Do not talk
Display Err - 2	Inflation error, pressure built up too	Check that air connector is correctly fitted,
	quickly or too slowly.	repeat measurement.
	Fault was recognised. Cuff is not	- Repeat measurement
	correctly fitted.	- Do not move arm
		- Do not talk
Display Err - 3	Deflation rate too high or too low.	Keep still during measurement. Call
		the Service Hotline if this fault occurs
		repeatedly.
Display 0	Start/Stop button accidentally	Switch unit off and then on again with
0	pressed while battery was being	Start/Stop button; if necessary.
	changed.	

Technical information **E** 

Failure encountered	Possible cause	Corrective action
\ <u>\\</u>	Irregular heartbeat, disruptive mo-	Repeat measurement after 3-5 minutes
-\\	vements, shaking, wobbling (objects),	rest. Possible factors (depending on the
	breathing in deeply, etc.	severity of the arrhythmia) affecting the
		measurement results must be discussed
		with your doctor.
Measured values are	Too little rest before the measure-	Repeat measurement after a break of
too high or too low	ment. Movement or talking during	approx. 3-5 mins. Observe instructions for
	measurement.	use P. EN-40.
	Has the proper cuff been selected?	Measure your arm circumference in the
	Cuffs that are too small result in	middle of the upper arm and compare this
	excessive blood pressure values.	with the details given on the cuff.
	Incorrect measurement position.	Repeat measurement, ensuring correct
		body posture. (Page EN-40)
Consecutive mea-	Blood pressure is not a fixed	
surements produce	quantity. It may vary in the upward	
different measured	or downward direction within a few	
values (readings).	minutes by more than 20 mmHg.	
Pulse values are too	Too little rest before the measure-	Repeat measurement after a break of
high or too low	ment. Movement or talking during	approx. 3-5 mins. Observe instructions for
	measurement.	use P. EN-40.
Measurement was	Unit detects a fault or that diastolic	Unit shows no fault. Repeat measurement.
interrupted and cuff	value is too low.	
first deflated and then inflated	Movement during measurement	If necessary, interrupt measurement, mea-
Lineir iiiitateu		sure again after a break of 5 minutes.

EN-53 EN-52

Technical informatio
La lecililicat illioilliatio

Failure encountered	Possible cause	Corrective action
No display after unit	Batteries inserted incorrectly.	Check position of batteries.
is switched on	Batteries flat.	Change batteries.
	Battery contacts dirty.	Clean battery contacts.
No display after	No measured values stored.	Take a new measurement.
Memory button is pressed		
Measurement	Batteries flat	Change batteries.
interrupted		
Date/time not shown	- Function accidentally deactivated	Activate date and time as described on P.
in display	- Function was not reactivated after	EN-42.
' '	battery replacement	

#### 2. Customer service

The unit may only be repaired by the manufacturer or by an expressly authorised organisation. Please contact:

UEBE Medical GmbH Phone: +49 (0) 9342/924040
Zum Ottersberg 9 Fax: +49 (0) 9342/924080
97877 Wertheim, Germany E-mail: info@uebe.com
Internet: www.uebe.com

Technical information

#### 3. Technical data

Unit type:	Digital automatic unit with electric pump for measuring blood pressure on upper arm	
Dimensions:	L = 100 mm x W = 152 mm x H = 61 mm	
Weight:	295 g not including batteries	
Display:	LCD display (liquid crystal display) for measured values and check displays	
Memory:	60 measured values (stored automatically) and average value (A)	
Measurement procedure:	Oscillometric determination of systole, diastole and pulse	
Reference procedure of clinical testing:	Auscultatory measurement	
Inflating pressure:	Initially up to approx. 180 mmHg, then in stages of 30 mmHg	
Pressure display range:	0-300 mmHg	
Measurement range	Systolic: 50-250 mmHg Diastolic: 40-150 mmHg Pulse: 40-160 pulses/min	
Error limits:	Blood pressure measurement: corresponds to EN 1060 Part 3 Pressure measurement: ± 3 mmHg Pulse measurement: ± 5%	
Serial number:	The unit bears a serial number which provides clear identification.	

Dower supply:	Battery type: 4 x 1.5 V AA alkali manganese (LR 6) or lithium (FR 6) round	
Power supply:	1 , ,,	
	cells, service life: more than 800 measurements in 2 years	
	Optionally stabilised type A1 mains adapter, output 6 VDC, minimum 600 mA	
Cuff:	Ring-cuff VWR (for arm circumference 22-43 cm)	
Operating conditions:	Ambient temperature 10 to 40 °C, rel. air humidity up to 85% not condensing, air pressure 700 to 1040 hPa	
Storage and transport conditions:	Ambient temperature -20 to 50 °C, rel. air humidity up to 85% not condensing	
Air deflation valve:	Electronically controlled	
Automatic switch-off	approx. 3 minutes after end of measurement	
IP code:	IP 20: Protected against solid foreign particles with a diameter of more than	
	12.5 mm, no protection against water.	
Expected service life:	5 years	

Internal power supply by battery

#### 4. Original spare parts and accessories:

The following original spare parts and accessories are available from specialist dealers:

 Ring-cuff 22-43 cm type VWR Part. no. 2505003 PZN 11312027

 Mains adapter type A1 Part. no. 2501020 PZN 00179482

Subject to technical modifications.

#### 5. Technical inspection / calibration check

Basically, a metrological inspection is recommended at intervals of 2 years. However, professional users in Germany are obligated to comply with the aforementioned according to "Regulation for Operators of Medical Devices".

This can be implemented either by UEBE Medical GmbH, an authority responsible for metrology, or authorised maintenance services. For this, please observe your national provisions.

Upon request, responsible authorities or authorised maintenance services receive a "Test instruction for metrological inspection" from the manufacturer.

Important: No modifications, e.g. opening the unit (except to replace the batteries) may be made to this unit without the manufacturer's prior permission.

EN-56

Classification:

Maintaining the unit

## **Technical information**

#### 6. Explanation of symbols

This product complies with the Council Directive 93/42/EC from 5 September 2007 regarding medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Degree of protection against electric shock: TYPE BF Please observe instructions for use

Relative humidity up to 85 %,

Serial number

Reference number = item number

Keep dry

Manufacturer

non-condensing

IP20

Protected against solid foreign particles with a diameter of more than 12.5 mm. no protection against water.

Equipment class II (double-insulated)

If using the optionally available mains adapter:


Internal fuse

Mains adapter specifications:

For indoor use only

Internal polarity positive

#### 7. Disposal

-20 °C-∕

Batteries and technical appliances must not be disposed of with domestic waste, but should be handed in at the appropriate collection and disposal points.

#### Maintaining the unit

- The unit contains sensitive parts and must be protected against strong variations in temperature, air humidity, dust and direct sunlight.
- Make sure when storing the unit that no heavy objects are left lying on the unit or on the cuff and that the air hose is not kinked.
- The unit is neither impact-resistant nor shock-proof. We recommend that you have the intactness and accuracy of the display checked after heavier falls or impacts.
- The unit is not waterproof.
- If the unit becomes dirty during use, use a soft, dry cloth to clean the unit. Do not use benzene, thinners or other strong solvents.
- Do not scrub or machine-wash the cuff. If the cuff gets dirty during use, use a synthetic cleaner and softly rub the surface.
- If the unit has been opened, it must be subjected to a metrological inspection conducted by an authorised organisation.
- Make sure that liquid cannot get into the air hose. Dry thoroughly.

## **G** Warranty

#### Warranty furnished

The device has been manufactured and tested with great care. However, in the unlikely event of a defect being detected after delivery, we provide warranty in accordance with the following terms and conditions:

During the warranty period of 2 years from the date of purchase we reserve the right either to repair any such defect at our expense or to supply a perfect replacement unit. The cost of returning the unit to our factory shall be borne by the sender. UEBE shall refuse to accept return deliveries that have not been paid for by the sender.

Excluded from the warranty are parts subject to normal wear and tear as well as damage caused by non-compliance with the instructions for use, improper handling (e.g. unsuitable power sources, breakages, leaking batteries) and/or disassembly of the unit by the purchaser. Furthermore, no claims for damages against us are substantiated by the warranty.

Warranty claims can only be advanced in the warranty period and by presenting proof of purchase. In the event of a warranty claim, the unit must be sent to the following address together with the proof of purchase and a description of the complaint:

UEBE Medical GmbH Service-Center Zum Ottersberg 9 97877 Wertheim Germany Warrai

englisl

EN-61

The cost of returning the unit to our factory shall be borne by the sender. UEBE shall refuse to accept return deliveries that have not been paid for by the sender.

The statutory claims and rights of the buyer against the seller (claims for defect, manufacturer's liability, for example) are not restricted by this warranty.

Please note:

In the event of a warranty claim it is essential to attach the proof of purchase.

EN-60

visocor® and UEBE® are internationally protected trademarks of

UEBE Medical GmbH Zum Ottersberg 9

97877 Wertheim

Germany

Phone: + 49 (0) 93 42 / 92 40 40 Fax: + 49 (0) 93 42 / 92 40 80

E-mail: info@uebe.com Internet: www.uebe.com

Subject to technical modifications.

Not to be reproduced, either in whole or in part.

© Copyright 2016 UEBE Medical GmbH