## INFECTIOUS MONONUCLEOSIS

### Intended use:

A rapid slide test for the qualitative and quantitative determination of infectious mononucleosis (IM) in human serum.

### Summary:

Paul and Bunell were the first to report that serum from a patient with infectious mononucleosis (IM) contained heterophile antibodies which agglutinated sheep erythrocytes. These heterophile antibodies react with an antigen which apparently is not responsible for their production. However, it was soon discovered that the test lacked specificity because the naturally occurring Forssman antibody found in serum from some individuals who apparently have not had recent infectious mononucleosis agglutinates unmodified sheep or horse erythrocytes. In 1937, Davidsohn employed a differential absorption procedure which removed the Forssman antibody yet retained the heterophile agglutination characteristic of IM. The Davidsohn modification added specificity but made the test time-consuming and cumbersome to perform. Therefore, the Davidsohn test has been relegated to the role of a reference method for diagnosis of IM. In the attempt to find a suitable alternative Bailey and Raffel discovered that bovine erythrocytes were more sensitive than sheep or horse erythrocytes for detecting IM heterophile antibodies. Since that time antigens, which have been extracted from bovine red cell membranes, have been used in various enzyme immunoassays (EIA) which are both highly sensitive and specific for heterophile antibodies associated IM.

### Test principle:

The Biocon IMslide Infectious Mononucleosis Latex Test provides a suspension of polystyrene latex particles which have been coated with partially purified glycoprotein from bovine red blood cells. The heterophile antibody associated with infectious mononucleosis binds to the corresponding antigenic determinants on the glycoprotein coated latex. This binding is evident by rapid agglutination of the latex. Due to the purification of the bovine red cell glycoprotein the coated latex is not agglutinated by Forssman or serum sickness antibodies at levels normally encountered therefore, no differential absorption is required.

### Reagent concentration:

#### R1:

Polystyrene Latex particles coated with partially purified glycoprotein from bovine red cells.

Preservative

#### R2/positive control:

Antigen containing solution

Preservative

#### R3/negative control:

Solution without Antigen

Preservative

### Preparation and stability:

The reagent and controls are stable up to the expiry date when stored at +2 to +8°C.

Do not freeze!

### Specimen:

Use fresh clear serum specimens.

Plasma, lipemic serum or microbial contamination may cause erroneous results. If the test cannot be performed immediately, store the specimen at +2 to +8°C for up to 24 hours. For longer storage, freeze the serum.

### Notes:

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Reagents containing sodium azide may combine with copper and lead plumbing to form highly explosive metal azides.

Dispose of reagent by flushing with large amounts of water to prevent azide buildup.

The positive and negative controls were prepared from human sera which have been tested using standard methods and found to be non-reactive for HbsAg and HIV antibodies. However, not test method can offer complete assurance that infectious agents are absent. Therefore all human specimens should be considered potentially infectious.

### Limitations – interference:

Although the Biocon IM slide Latex Reagent is highly sensitive and specific, a diagnosis of infectious mononucleosis should not be made on the basis of a positive test result without the support of patient history and hematological or other clinical evidence. Similarly, a negative test result cannot completely rule out infectious mononucleosis. Incubation of the test for longer than the recommended time or microbial contamination may cause false positive reactions. Apparent false positive reactions have been associated with sera from patients with other disease such as rheumatoid arthritis, certain respiratory infections, leukaemia, Burkitt's lymphoma and serum sickness. Although most patients develop heterophile antibodies within 3 weeks of the onset of symptoms, occasional patients may take several months to develop detectable levels. If the Biocon IMslide IM Latex Test is negative in the presence of strong evidence suggesting a diagnosis of infectious mononucleosis, repeat testing on samples obtained at intervals of several days will generally reveal development of the heterophile agglutinin. Some patients with hematological and clinical evidence of infectious mononucleosis remain persistently negative. A single heterophile antibody titer cannot be interpreted as an indication of the stage or severity of the disease. However, titrations on sequential may be useful in following the course of the disease in an individual patient.

### Testing procedure:


#### Qualitative test:

1. Bring reagents and specimens to room temperature before use.
2. Place one drop (1 x 50 µl) of the positive control on field # 1 of the reaction slide.
3. Place one drop (1 x 50 µl) of the negative control on field # 2. Using a pipette, place one drop (50 µl) of each test specimen on successive fields.
4. Gently resuspend the Latex Reagent and add one drop (35 µl) to each test field. Use stick to spread reaction mixture over entire test field.
5. Rotate the slide for 2 minutes and read **immediately** under direct light.

#### Quantitative test:

1. Bring reagents and specimens to room temperature before use.
2. Using saline, dilute the specimens 1:2, 1:4, 1:8, 1:16 or as needed.
3. Place one drop (50 µl) of each dilution on successive fields of the reaction slide. Gently resuspend the IM Latex Reagent and add one drop (50 µl) to each test field. Use stick to spread reaction mixture over entire test field.
4. Rotate the slide for 2 minutes and read **immediately** under direct light.

### Interpretation:

### Results:

#### Qualitative test:

A negative reaction is indicated by a uniform milky suspension with no agglutination as observed with the negative control (2). A positive reaction is indicated by any observable agglutination in the reaction mixture. The specimen reaction should be compared with the positive control (1).

#### Quantitative test:

The titer of the serum is the reciprocal of the highest dilution which exhibits a positive reaction.

# IMslide

## INFECTIOUS MONONUCLEOSIS

### **Sensitivity and method comparison:**

Serum specimens from 285 individuals that had been submitted to clinical laboratories by physicians for IM testing were examined. The Biocon IMslide and another commercial Red Cell IM Test Kit were used to evaluate the specimens. One hundred thirty-two (132) specimens were found positive by both assays. The remaining 153 specimens gave negative results using both products. These data indicate that both sensitivity and specificity of the Biocon IMslide Latex are 100%.

### **Precision:**

In a study on precision, a panel of 10 serum samples with IM heterophile antibody titers from 1 to 256 were tested on 10 consecutive days by the Quantitative Method (100 determinations). No determinations gave more than a 2-fold difference from the mean titer for a sample.

### **Quality control:**

Biocon IMslide **R2/positive control** and **R3/negative control** should be included in each test series. The **R2** should produce strong agglutination and the **R3** should yield a smooth suspension with no agglutination.

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