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## Introduction

The Wobblesmart® is used to improve balance, posture, and coordination. Training on the unstable platform requires balance and a strong posture. It is mainly used in prevention (children and seniors), in rehabilitation after an injury, and in sports in order to enhance performance.

Studies prove that targeted balance training can reduce the frequency of falls and secondary injuries. Due to higher demands placed on supporting muscles, the balance trainer can also improve your performance with conventional exercises.

## Wobblesmart®

### Materials and characteristics

The Wobblesmart® consists of wood and offers maximum safety for balance training due to its plastic foot and non-slipping surface area. Even people with bigger feet can stand two-legged on the 39 cm diameter platform and find a stable stance. The height respectively difficulty can be adjusted between 65 mm (level 1) and 80 mm (level 6).

Level	Instability
1	low
2	moderate
3	strong
4	very strong
5	special strong
6	extremely strong

### **Precautions**

Every time you train, first verify that the plastic foot of the Wobblesmart® is resting firmly on the surface below it and cannot slide. Avoid training on slippery floors or rugs. The Wobblesmart® may be placed onto a gymnastics mat in order to prevent slipping. However, this does decrease the level of difficulty for exercises. Persons suffering from pronounced balance, posture, or sensory disturbances and persons with high risk of falls (caused by medications, osteoporosis, or reduced muscle strength) should only train under the instruction of a physician, physiotherapist, or specialised trainer.

- » Remove sharp or pointed objects from the immediate vicinity of the training area.
- » Regular cleaning with warm water and soap is recommended. For disinfectant use the ARTZT disinfectant spray.
- » The training can either be performed with shoes, bare-foot or with socks. We recommend the usage without shoes.

## Training tips

For effective training to improve and enhance your balance, progress from basic to advanced exercises:

- » Use a sturdy, fixed object for support while performing the exercises (for example, you can hold on to a wall or a sturdy chair).
- » Train close to the chosen support object so that you can hold yourself up if necessary.
- » Progress to train without any support at all.

Start the training with Wobblesmart® on level 1, so that the plastic foot offers the maximum contact face to the floor whilst the instability is low.

Not before feeling absolutely confident and being able to hold the balance during a complete exercise sequence, you may increase the level of difficulty by rotating the foot clockwise.

The increase in training intensity should be tailored to individual performance so that the exercises are always performed safely, correctly, and with good posture and coordination.

### Adjustable foot (Level 1-6)

Health-oriented, preventive training focuses on the improvement of the general sense of balance. Beginners should select 6 to 8 of the following exercises in order to achieve this goal. Static balance exercises (exercises without movement) should be held 3 times for 10-30 seconds at a time. Dynamic exercises (motion sequences, stepping, etc.) should be repeated 10-15 times in 1-3 sets. You should never experience EXTREMELY heavy strain while training. Once you have established a foundation and the exercises come easily, you can increase the level of difficulty, number of repetitions (sets), and training frequency. Please observe the following principles:

- » Maintain correct posture and conduct movements in a slow, controlled manner in the initial position, during the exercise, and at the end.
- » Training should be pain-free. Take a break if an exercise becomes painful. You should consult a healthcare provider if pain occurs repeatedly.
- » Take short breaks or change sides between static exercises. You may also alternate between dynamic and static exercises.
- » Examine each exercise in detail and visualise the process.
- » ALWAYS exercise both the right and left side of the body even when the illustration only shows one side.

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Stand on the Wobblesmart® and find a stable position (without wobbling).

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Rotate your head slowly rightwards and then leftwards.

**Zwei-Bein-Stand | Two-legged stance**  
**Augen zu | Closed eyes**

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Find a stable position and close your eyes.

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Stand on the Wobblemart® and find a stable position.

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Raise your right leg to the front and return to the initial position. Alternate leg.

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Shift your weight to the right leg and lift the left heel backwards. Find a stable position. Alternate leg.

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Find a stable stance. Keep your gaze directed ahead.

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Bend your knees to a squat while moving your extended arms to the front. Now slowly return to the initial position.

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Stand on the Wobblesmart® and find a stable stance.

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Alternately step down with the left and right foot, with only your toes touching the ground. Always resume the centred position in between steps.

**Ein-Bein-Stand | One-legged stance**  
**Bein- und Armbewegungen | Leg and arm movements**

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Stance behind the Wobblesmart®. Step carefully with one foot on the Wobblesmart® and find a stable position. Alternate legs. If you are able to complete this exercise easily you can lift the arms additionally forward and up.

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Step carefully on the Wobblesmart® with one foot and find a stable position. Move the free leg with flexed knee to the side. Move your arms back slightly from their forward position. Alternate leg.

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Stance on the side next to the Wobblesmart®. Step carefully with one foot on the Wobblesmart® and lift the other foot from the ground. Find a stable position.

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Just like A, but now move the free foot to the knee of the opposite leg. The arms are slightly bent. Keep the position. Alternate leg.

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While standing one-legged on the Wobblesmart®, lift an object (e.g. a towel or ARTZT vitality® Latexfree exercise band) from the ground and then put it back.

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Get into the push-up position whereas the knees work as points of support.

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Perform a push-up by bending and extending your elbows and keeping them close to your trunk, while trying to stabilize the Wobblesmart®.

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Get into the push-up position with extended legs. Perform the push-ups.

# UPGRADE YOUR VITALITY!

**ARTZT**  
vitality

## Produktprogramm ARTZT vitality® | Product range ARTZT vitality® Medical Fitness - Made in Germany

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- **Flowin® Personal Training**
- **Gymstick™ Health**
- **BOSU® Balance Training**

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ARTZT vitality® is a registered trademark of Ludwig Artzt GmbH. Unauthorized use is strictly prohibited.  
© 2009 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.

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© 2009 Ludwig Artzt GmbH. Alle Rechte geschützt. Gedruckt in Deutschland.