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Introduction

As specialist for fitness and health, Ludwig Artzt GmbH has been building a bridge between sports science and the fitness world for over 25 years now; the product range especially focuses on functional and flexible small devices. The brand ARTZT vitality® offers professionals and home users a complete range of high-quality products and modern, efficient training concepts for professional and recreational athletes.

Permanent, close dialogue with experts is the mainstay of Ludwig Artzt GmbH's success. Our in-house (advanced) training institute is known as one of the best qualified seminars on health and fitness issues; it particularly focuses on conveying knowledge based on the findings of sports science.

With the Massage Roll by ARTZT vitality®, an effective and mobile training device that stimulates active regeneration is being added to the product range. The ARTZT vitality Massage Roll mobilises, repairs the muscle strain and results the muscle to relax. It also promotes coordination through a "special" roller training.

For sports or everyday use – the ARTZT vitality® Massage Roll offers active compensation for the cumulative effects of inappropriate and excessive strain that are typical for our modern, high-speed society.

Massage Roll

Materials and characteristics

The ARTZT vitality® Massage Roll offers flexible training for anyone. It features high-quality workmanship with the following positive characteristics:

Mobile & lightweight

with a weight of just 145 g and compact dimensions (30 cm long; 15 cm; hole diameter: 5 cm), the Massage Roll can be used for professional training just about anywhere (indoor and outdoor).

Straightforward & versatile

Thanks to easy handling, the Massage Roll is suitable for everyday use, therapy and sports.

Efficient

For self myofascial release in particular, the user can achieve a normal muscle tension with the help of his or her own body weight – without the assistance of a therapist.

Robust & abrasion-resistant

Easy to clean, can be sterilised, inherent thermal stability up to 110 degrees Celsius.

Top quality & safe

The recyclable material (EPP) is free of propellant gases and other chemical propellants, odourless and not soluble in water. As a result, using the Massage Roll is guaranteed to be physiologically harmless.

Made in Germany

Possible Applications

The ARTZT vitality® Massage Roll is much more versatile than it seems at the first sight. It is suitable for use in physiotherapy, competitive sports, amateur sports and fitness sports. The Roll has many uses:

- » Mobilising and relaxing muscle adhesions with the help of myofascial release techniques (= loosening / massaging agglutinated muscle fibre bundles, accompanied by improved functionality of the musculature)
- » Active regeneration with improved circulation
- » Ideal warm up before sports and regeneration after
- » Reduction of muscle tension and stiffness
- » Supportive cellulite treatment / lymph drainage
- » Coordinative strength training, especially as an additional part of Pilates training*
- » Balance exercises


These instructions focus on describing the myofascial self-relaxation technique. The difference between this and other techniques is that the practitioner utilises his or her own body weight in order to complete the exercises, rather than having to rely on the help of another person.

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Active Principle
Golgi Tendon Organ

Myofascial self-relaxation with the Massage Roll is made possible by the principle of autogenous inhibition. Here the Golgi Tendon Organ plays a major role since it responds to changes in muscle tension. When tension increases, the Golgi Tendon Organ responds and relaxes it. The Massage Roll is an ideal tool for stimulating the Golgi Tendon Organ while simultaneously achieving a reduction in muscle tension.

Abnormal fascia changes can lead to chronic pain, reduced flexibility and a decrease in performance. Myofascial self-relaxation with the Massage Roll constitutes an effective, affordable, straightforward method to treat myofascial scar tissue and adhesion. The Massage Roll is the ideal regeneration device for self-massage since it can help you to release fascial adhesion and break down scar tissue.

Safety information

The durability of and safe training with the Massage Roll are largely dependent on proper handling and correct use:

- » Remove sharp-edged or pointy objects from the training area and ensure you have enough space, since some of the movements are expansive.
 - » If the Massage Roll gets dirty after intensive use by various people, it can be cleaned easily using a damp cotton cloth. Soap and water may be used for cleaning very dirty rolls. Only use the roll after it has dried completely.
 - » Avoid training on slippery or rough floors. The Massage Roll should roll and not slight.
 - » The Massage Roll is designed for fitness, sports and therapy applications. It may not be used for other purposes.
 - » During coordinative balance exercises in a standing position, the point load on the Massage Roll must not exceed 100 kg.
- » Improper or careless use of the Massage Roll may lead to injuries. As a result, only safe and pain-free exercises should be completed.
 - » The Massage Roll is not intended as a toy for children! Children must be supervised during use!
 - » If you suffer from osteoporosis, intervertebral disc lesions, fibromyalgia, thrombosis or elevated blood pressure, talk to your doctor or therapist before using this product.

Training frequency

An active self-massage can be completed several times a day. Each training unit should be conducted for approximately 2 minutes at 30-second intervals. 30 seconds of active use of the Roll alternate with pauses of 30 seconds.

When using the Massage Roll for coordinative strength exercises, the frequency recommendations for health-oriented strength training apply (1-3 sets with 15-20 repetitions each). The exercises described below focus on active self-massage.

Training Preparation and Execution

Please look over the following principles:

- » Start your training slowly and do not put excessive demands on yourself.
 - » Ensure that you have sufficient space for training.
 - » Maintain correct posture and move in a slow controlled way in the initial position, during the exercise, and at the end.
- » Roll the Massage Roll over the tense muscle group. When you reach a tense, painful point, pause for a few seconds and then continue rolling. By supporting the arms or the inactive leg (during exercises while lying down), you can vary the pressure applied to the Roll while simultaneously increasing torso stability. During exercises conducted while standing up, you can regulate the massage effect yourself at any time through the pressure of the upper body, depending on your sensitivity.
 - » Pain may result the first one or two times when the Massage Roll is used for the treatment of muscle problems (comparable to a massage with pronounced muscle tension). However, pain should decrease and then disappear rapidly with regular use. If pain recurs over an extended period of time, please consult a doctor.
 - » The exercises with the Massage Roll are illustrated with the initial and final position, respectively. Study each exercise in detail and visualise the process.

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Stand with your back to the wall and hold the Massage Roll between the wall and lower back. The feet are approximately 40 cm away from the wall.

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Now bend your legs and do a squat without reducing pressure on the Massage Roll.

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Stand with your back to the wall and hold the Massage Roll between the wall and back at the level of the thoracic spine (chest height). The feet are approximately 40 cm away from the wall

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Stand with your side to the wall and hold the Massage Roll between the wall and waist. The feet are approximately 20 cm away from the wall, the knees are slightly bent.

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Stand with your side to the wall and hold the Massage Roll between the wall and upper arm. The feet are approximately 20 cm away from the wall, the knees are slightly bent.

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Stand with your side to the wall and hold the Massage Roll between the wall and forearm, shortly above the wrist. The feet are approximately 20 cm away from the wall, the knees are slightly bent.

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Lie on the floor, on your side. Support yourself on your left elbow (the upper body is not touching the floor) and place the Massage Roll below your left thigh (shortly above the knee). The right foot and the right hand are supported in front of the body.

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Now push the upper body towards your foot with the help of your arms. The upper body remains in an upright, sideways position during the entire movement.

Oberschenkelrückseite / Ischiokrurale Muskulatur | Back of the thigh / ishiocrural musculature

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Sit on the floor with your legs extended. Support your arms behind the body and place the Massage Roll underneath your thighs (shortly below the buttocks). Your feet are suspended above the floor.

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Now pull your body towards your hands with the help of your arms. The upper body remains upright during the entire movement.

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
Sit on the floor with your legs extended. Support your arms behind your body and place the Massage Roll underneath the calves (shortly above the ankles). Your feet are suspended above the floor.

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Now push the body towards your feet with the help of your arms. The upper body remains upright during the entire movement.

UPGRADE YOUR VITALITY!

ARTZT
vitality

Produktprogramm ARTZT vitality® | Product range ARTZT vitality® Medical Fitness - Made in Germany

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- **Flowin® Personal Training**
- **Gymstick™ Health**
- **BOSU® Balance Training**

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ARTZT vitality® is a registered trademark of Ludwig Artzt GmbH. Unauthorized use is strictly prohibited.
© 2009 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.

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© 2009 Ludwig Artzt GmbH. Alle Rechte geschützt. Gedruckt in Deutschland.