


Chek-Stix UA Diptube Liquid Quality Control

[siemens-healthineers.com/chek-stix](https://www.siemens-healthineers.com/chek-stix)



Quality control made effortless

As a leader in urinalysis, Siemens Healthineers strives to deliver quality solutions you can trust. Healthcare and laboratory professionals rely on Siemens Healthineers to ensure consistent, accurate results and help streamline operations. Building on that legacy of leadership and trust, Siemens Healthineers is proud to deliver the latest innovation in its urinalysis portfolio—introducing Chek-Stix® UA Diptube Liquid QC solution.



Key benefits:



Liquid-ready solution means less time prepping quality control solution and more time running patient samples.



Barcode-enabled for quick, accurate data input.



Ability to perform multiple dips (20 times) in the same tube.



30-day room temperature stability.



Supports our broad urine chemistry menu, including albumin, creatinine, blood, glucose, and hCG.



Validated for use with all Siemens Healthineers urinalysis devices and reagents.



22 months' shelf life from manufacture date.


Chek-Stix UA Diptube Liquid QC

Streamline your daily operations with a liquid-ready urinalysis solution for use with Siemens Healthineers point-of-care urinalysis instruments and reagents, including 12 parameters and hCG.



Product Configuration

Chek-Stix UA Diptube Liquid QC contains both Level 1 (negative for all chemistry parameters and positive for hCG) and Level 2 (positive for all chemistry parameters and negative for hCG), packaged in one box. Both Level 1 and Level 2 diptube controls are ready-to-use liquid quality controls for use with the CLINITEK Status®+ Analyzer, CLINITEK Advantus® Analyzer, and visual-read urine strips.



Analyte Profile

pH
Protein
Glucose
Ketone
Bilirubin
Blood (Hemoglobin)
Nitrite
Urobilinogen
Specific Gravity
Leukocyte Esterase
Creatinine
Albumin
A:C (Albumin-to-Creatinine Ratio)
P:C (Protein-to-Creatinine Ratio)
hCG

Ordering Information

SMN	Each kit contains
11561694	Level 1 3x15mL Level 2 3x15mL

Chek-Stix, CLINITEK Advantus, CLINITEK Status, and all associated marks are trademarks owned by or licensed to Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
Phone: +49 9191 18-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics Inc.
Point of Care Diagnostics
2 Edgewater Drive
Norwood, MA 02062-4637
USA
Phone: +1 781-269-3000